# Evaluation of the Performance of CPAP Heated Humidifiers for use in Sleep Apnea Therapy; A Comparative Study of Humidification Effectiveness

R. Virag, *Trifid Medical Group, United States*. In submission

### **KEY FINDINGS:**

- ThermoSmart<sup>®</sup> is the only CPAP humidifier in this study to deliver average absolute humidity greater than 27 mg/L.
- **ThermoSmart**<sup>®</sup> is the only CPAP humidifier in this study to deliver high levels of humidity, without causing any condensate.

#### AIM:

To compare the temperature and humidity performance of different CPAP humidifiers under simulated real-use conditions.

#### **METHODS:**

- Four humidified CPAP devices from different manufacturers were side-by-side bench-tested over seven hours.
- Test conditions included a CPAP pressure of 12 cmH<sub>2</sub>O, 40 L/min airflow, maximum humidifier settings, ambient temperature of 67.3 °F (19.6 °C), and ambient humidity of 9.3 mg/L.
- End-of-hose temperatures were recorded hourly for each device.
- Two methods were used to record humidity: a) Gravimetric method using ISO 8185; b) Hygrometry.



## **KEY:**

185045104 REV B en

FPH – Fisher & Paykel Healthcare SleepStyle<sup>™</sup>608 with ThermoSmart<sup>®</sup> Respironics – REMstar<sup>®</sup> Pro M-Series with heated humidifier REMstar<sup>®</sup> is a registered trademark of Respironics Inc **ResMed 1** – S8 Elite<sup>™</sup> HumidAire 3i<sup>™</sup> version 1 **ResMed 2** – S8 Elite<sup>™</sup> HumidAire 3i<sup>™</sup> version 2 S8 Elite<sup>™</sup> and HumidAire 3i<sup>™</sup> are trademarks of ResMed Ltd

#### CONCLUSION:

• ThermoSmart<sup>®</sup> outperformed the other devices in this study in all areas of heat and humidity delivery, with zero condensation.



