

Table of Contents

Intended Use	2
Warnings	2
Cautions	3
Contraindications	4
Symbol Key	4
System Contents	4
System Overview	5
Control Buttons	6
Installing the Air Filters	6
Connecting the Breathing Circuit	7
Where to Place the Device	7
Supplying AC Power to the Device	7
Setup the Integrated Humidifier	8
Display	9
Starting the Device	10
Heated Humidifier Feature	11
Ramp Feature	11
C-Flex Feature	12
Navigating the Display Screens	12
Viewing Data on the Patient Data Screens	12
System Error Screen	13
Device Alerts	14
Troubleshooting	14
Accessories	15
Adding Supplemental Oxygen	15
Traveling with the System	15
Cleaning the Device	15
Cleaning or Replacing the Filters	15
Cleaning the Tubing	16
Cleaning the Water Chamber	16
Service	16
Specifications	17
Disposal	18
How to Contact Respironics	18
EMC Information	19
Limited Warranty	Back Page

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician.

Read all instructions before using the SleepEasy with C-Flex device with integrated humidifier.

Intended Use

The Respironics SleepEasy with C-Flex device delivers continuous positive airway pressure (CPAP) therapy for the treatment of Obstructive Sleep Apnea in spontaneously breathing patients weighing >66 lbs. It is intended for use in the home or hospital environment.

IMPORTANT: Your home care provider will make the correct pressure settings according to your health care professional's prescription. Several accessories are available to make your OSA treatment with the SleepEasy system as convenient and comfortable as possible. To ensure that you receive the safe, effective therapy prescribed for you, use only Respironics accessories.

Warnings

A warning indicates the possibility of injury to the user or the operator.

- This manual serves as a reference. The instructions in this manual are not intended to supersede the health care professional's
 instructions regarding the use of the device.
- The operator should read and understand this entire manual before using the device.
- This device is not intended for life support.
- The device should be used only with masks and connectors recommended by Respironics or with those recommended by the health care professional or respiratory therapist. A mask should not be used unless the device is turned on and operating properly. The exhalation port(s) associated with the mask should never be blocked. Explanation of the Warning: The device is intended to be used with special masks or connectors that have exhalation ports to allow continuous flow of air out of the mask. When the device is turned on and functioning properly, new air from the device flushes the exhaled air out through the mask exhalation port. However, when the device is not operating, not enough fresh air will be provided through the mask, and exhaled air may be rebreathed. Rebreathing of exhaled air for longer than several minutes can in some circumstances lead to suffocation.
- Do not use this device without the water chamber installed.
- This device contains small parts which could present a choking hazard.
- · Do not use extension cords with this device.
- If oxygen is used with the device, turn on the device before you turn on the oxygen flow. Turn the oxygen off before turning the
 device off. Also, the oxygen flow must be turned off when the device is not in use. Explanation of the Warning: When the
 device is not in operation and the oxygen flow is left on, oxygen delivered into the tubing may accumulate within the device's
 enclosure. Oxygen accumulated in the device enclosure will create a risk of fire.
- When using oxygen with this system, a Respironics Pressure Valve must be placed in-line with the patient circuit after the oxygen source. Failure to use the pressure valve could result in a fire hazard.
- Do not use the device in the presence of a flammable anaesthetic mixture in combination with oxygen or air, or in the presence
 of nitrous oxide.
- When using oxygen with this system, the oxygen supply must comply with local regulations for medical oxygen.
- Oxygen accelerates fires. Keep the device and the oxygen container away from heat, open flames, any oily substance, or other sources of ignition. Do not smoke in the area near the device or the oxygen.
- If administering fixed-flow supplemental oxygen, the oxygen concentration may not be constant. The inspired oxygen
 concentration will vary, depending on the CPAP setting, patient breathing pattern, and leak rate. Substantial leaks around the
 mask may reduce the inspired oxygen concentration to less than the expected concentrations. Appropriate patient monitoring
 should be implemented.
- · Do not connect the device to an unregulated or high pressure oxygen source.
- · Do not use the device near a source of toxic or harmful vapors.
- Do not use this device if the room temperature is warmer than 95° F. If the device is used at room temperatures warmer than 95° F, the temperature of the airflow may exceed 106° F. This could cause irritation or injury to your airway.
- Do not operate the device in direct sunlight or near a heating appliance because these conditions can increase the temperature of the air coming out of the device.
- Contact your health care professional if symptoms of sleep apnea recur.
- If you notice any unexplained changes in the performance of this device, if it is making unusual or harsh sounds, if it has been
 dropped or mishandled, if the enclosure is broken, or if water has entered the enclosure, disconnect the power cord and
 discontinue use. Contact your home care provider.
- Repairs and adjustments must be performed by Respironics-authorized service personnel only. Unauthorized service could cause injury, invalidate the warranty, or result in costly damage.
- Inspect the power cord often for any signs of damage. Discontinue use and replace if damaged.

- · To avoid electric shock, unplug the device before cleaning it. DO NOT immerse the device in any fluids.
- Pins of connectors identified with the ESD warning symbol () should not be touched. Connections should not be made to these connectors unless ESD precautionary procedures are used. Precautionary procedures include methods to prevent build-up of electrostatic charge (e.g., air conditioning, humidification, conductive floor coverings, non-synthetic clothing), discharging one's body to the frame of the equipment or system or to earth or a large metal object, and bonding oneself by means of a wrist strap to the equipment or system or to earth.
- Never operate the device if any of the parts are damaged, or if it is not working properly. Do not use the device if the water chamber is leaking or damaged in any way. Have any damaged parts replaced before continuing use.
- · Never touch the heater plate unless the unit is unplugged and the plate has cooled down.
- · Allow the water in the chamber to cool to room temperature before removing the chamber from the humidifier.
- While the device is in operation, the soft flap on the flapper valve must move freely and close off the opening on top of the water chamber. Replace the flapper valve if the flap is damaged or not intact.
- · When installing the water chamber onto the system, do not allow any water to spill into the device.
- Empty and clean the water chamber daily to prevent mold and bacteria growth.
- · When installing the water chamber onto the system, make sure the soft valve is first attached to the water chamber.
- If the device is used by multiple persons (e.g., rental devices), a low-resistance, main flow bacteria filter should be installed in-line between the device and the circuit tubing to prevent contamination.
- The exhalation device (e.g., Whisper Swivel II) or exhalation port (on masks with an integrated exhalation port) is designed to
 exhaust CO, from the patient circuit. Do not block or seal the ports on the exhalation device.
- If you are using a full face mask (a mask covering both your mouth and your nose), the mask must be equipped with a safety (entrainment) valve. You must ensure that the entrainment valve is functioning properly.
- · Route the wires to avoid tripping.
- This device is activated when the power cord is connected. Pressing the START/STOP button (1) turns the airflow and humidifier (if activated) on or off.

Cautions

A Caution indicates the possibility of damage to the device.

- The device may only be operated at temperatures between 41° F and 95° F.
- If this device has been exposed to either very hot or very cold temperatures, allow it to adjust to room temperature (operating temperature) before starting therapy. Do not operate the device outside of the operating temperature range shown in the Specifications section.
- · Do not immerse the device or allow any liquid to enter the enclosure, the inlet filter, or any opening.
- Do not place the device in or on any container that can collect or hold water.
- Make sure that the drain holes on the bottom of the device are not blocked.
- Condensation may damage the device. Always allow the device to reach room temperature before use.
- · A properly installed, undamaged reusable foam inlet filter is required for proper operation.
- Tobacco smoke may cause tar build-up within the device, which may result in the device malfunctioning.
- If fluids are spilled onto the heater plate, unplug the power cord from the AC wall outlet and allow the device to drain and dry before using.
- Take precautions to protect furniture from water damage.
- · Use only distilled water in the chamber. The use of additives or chemicals in the chamber is prohibited.
- Do not fill water chamber above the fill line indicated on the side of the water chamber. Damage to the therapy device may
 occur.
- Make sure the protective cover is securely fastened on top of the water chamber before operating the device.
- Do not place the device directly onto carpet, fabric, or other flammable materials.
- Do not allow the water chamber to sit for any length of time after it has been filled with water. Immediately install the chamber
 in the device. Allowing water to sit in the chamber (when the chamber is not installed in the device) may cause the chamber to
 separate from the bottom plate and may cause water leakage.
- · Avoid moving the device when the water chamber has water in it.
- · Operating the device with a dirty filter may keep the system from working properly and may damage the device.
- Dirty inlet filters may cause high operating temperatures that may affect device performance. Regularly examine the inlet filters as needed for integrity and cleanliness.
- · Never install a wet filter into the device. You must ensure sufficient drying time for the cleaned filter.

Contraindications

When assessing the relative risks and benefits of using this equipment, the clinician should understand that this device can deliver pressures up to 20 cm H_2O . In the event of certain fault conditions, a maximum pressure of 30 cm H_2O is possible. Studies have shown that the following pre-existing conditions may contraindicate the use of CPAP therapy for some patients:

- · Bullous Lung Disease
- · Pathologically Low Blood Pressure
- Bypassed Upper Airway
- Pneumothorax
- Pneumocephalus has been reported in a patient using nasal Continuous Positive Airway Pressure. Caution should be used
 when prescribing CPAP for susceptible patients such as those with: cerebral spinal fluid (CSF) leaks, abnormalities of the
 cribriform plate, prior history of head trauma, and/or pneumocephalus. (Chest 1989; 96:1425-1426)

The use of positive airway pressure therapy may be temporarily contraindicated if you exhibit signs of a sinus or middle ear infection. Not for use with patients whose upper airways are bypassed. Contact your physician if you have any questions concerning your therapy.

Symbol Key

The following symbols may appear on the device:

Symbol	Definition	Symbol	Definition
[]i	Consult accompanying instructions for use.		Electrostatic Discharge
RONLY	Caution! U.S. federal law restricts this device to sale by or on the order of a physician.		Hot Surface
沈	Type BF Applied Part	\triangle	Caution
	Class II (Double Insulated)	10101	Serial Interface
IPX1	Drip Proof Equipment	7	Maximum Fill Line

System Contents

Your SleepEasy system includes the following items:

- · SleepEasy device with integrated humidifier
- AC power cord
- · Flexible tubing
- · Protective cover
- Funnel
- Reusable gray foam filter and disposable ultra-fine filter
- · Travel bag

Note: If any of these items are missing, contact your home care provider.


System Overview

The SleepEasy device is a sleep apnea system that delivers Continuous Positive Airway Pressure (CPAP). CPAP maintains a constant level of pressure throughout the breathing cycle.

When prescribed for you, the device provides several special features to help make your therapy more comfortable. The ramp function allows you to lower the pressure when you are trying to fall asleep. The air pressure will gradually increase until your prescription pressure is reached. You also have the option of not using the ramp feature at all.

Additionally, the C-Flex comfort features provide you with pressure relief when you exhale during therapy.


Several accessories are also available for use with your SleepEasy device. Contact your home care provider to purchase any accessories not included with your system. The following diagram illustrates many of the device features.



The features are described in the following table.

Device Feature	Description
Heater Plate	Warms the water in the water chamber.
Air Outlet Port (conical, 22 mm)	Connect the flexible tubing here.
Water Chamber	The removable water chamber holds the water for humidification.
Fill Line	This indicates the maximum water level for safe operation.
Protective Cover	Clip on top of the water chamber to prevent debris from entering the water chamber.
Flapper Valve	This valve helps prevent water from splashing into the device.
Power Inlet	Connect the power cord here.
Medical Equipment Note (bottom of unit)	For ease at airport security stations, there is a note on the bottom of the device stating that it is medical equipment. It may help if you also take this manual with you when you travel.
Filter Area (back of unit)	A reusable, gray foam filter must be placed in the filter area to screen out normal household dust and pollens. An optional, white ultra-fine filter can also be used for more complete filtration of very fine particles.

Control Buttons



These buttons are described below.

Виттом	Description
0	Starts the airflow and places the device into Active state, or stops the airflow, and places the device into Standby state. It is also used to exit any display screen.
XX.	Controls the humidifier functions.
	When the airflow is on, this button allows you to activate or restart the ramp function. Ramp lowers the airflow pressure and then gradually increases it, allowing you to fall asleep more easily.
Flex	Controls the C-Flex functions. This allows you to adjust the level of air pressure relief that you feel when you exhale during therapy.
+	Performs next screen navigation or increases setting.
-	Performs previous screen navigation or decreases setting.

The START/STOP, HEAT, RAMP and FLEX buttons are lit by LED backlights. The START/STOP LED will always be on when power is applied to the device. The HEAT LED is on when the humidifier heater plate is active, and off when it is not. The RAMP LED is on when the ramp function is active and off when it is not. The HEAT, RAMP and FLEX LEDs will flash when their respective settings are being adjusted.

Installing the Air Filters

CAUTION: A properly installed, undamaged reusable gray foam filter is required for proper operation.

The device uses a gray foam filter that is washable and reusable, and an optional white ultra-fine filter that is disposable. The reusable filter screens out normal household dust and pollens, while the optional ultra-fine filter provides more complete filtration of very fine particles. The gray reusable filter must be in place at all times when the device is operating. The ultra-fine filter is recommended for people who are sensitive to tobacco smoke or other small particles.

A reusable gray foam filter and a disposable ultra-fine filter are supplied with the device. If your filters are not already installed when you receive your device, you must at least install the reusable gray foam filter before using the device.

To install the filter(s):

- 1. If you are using the white disposable ultra-fine filter, insert it into the filter area first, mesh-side facing in, towards the device.
- 2. Insert the gray foam filter into the filter area after the ultra-fine filter.

Note: If you are not using the white disposable filter, simply insert the gray foam filter into the filter area.

Connecting the Breathing Circuit

To use the system, you will need the following accessories in order to assemble the recommended circuit:

- Respironics interface (nasal mask or full face mask) with integrated exhalation port, or Respironics interface with a separate exhalation device (such as the Whisper Swivel II)
- · Respironics 6 ft. flexible tubing
- Respironics headgear (for the mask)

WARNING: If the device is used by multiple persons (e.g., rental devices), a low-resistance, main flow bacteria filter should be installed in-line between the device and the circuit tubing to prevent contamination.

To connect your breathing circuit to the device, complete the following steps:

1. Connect the flexible tubing to the air outlet on the top of the device.

Note: If required, connect a bacteria filter to the device air outlet, and then connect the flexible tubing to the outlet of the bacteria filter.

Note: Do not use the heated humidifier feature when using a bacteria filter. No water should be in the humidifier chamber when using the bacteria filter.

Note: The bacteria filter is recommended to protect the patient, care provider and equipment from the transference of a virus or bacteria through the breathing circuit.

2. Connect the tubing to the mask. Refer to the instructions that came with the mask.

WARNING: The exhalation device (e.g., Whisper Swivel II) or exhalation port (on masks with an integrated exhalation port) is designed to exhaust CO₂ from the patient circuit. Do not block or seal the ports on the exhalation device.

WARNING: If you are using a full face mask (a mask covering both your mouth and your nose), the mask must be equipped with a safety (entrainment) valve. You must ensure that the entrainment valve is functioning properly.

3. Attach the headgear to the mask if necessary. Refer to the instructions that came with your headgear.

Where to Place the Device

Place the device on a firm, flat surface somewhere within easy reach of where you will use it at a level lower than your sleeping position. Make sure the filter area on the back of the device is not blocked by bedding, curtains, or other items. Air must flow freely around the device for the system to work properly. Make sure the device is away from any heating or cooling equipment (forced air vents, radiators, air conditioners).

CAUTION: Do not place the device directly onto carpet, fabric, or other flammable materials.

CAUTION: Do not place the device in or on any container that can collect or hold water.

CAUTION: Make sure that the drain holes on the bottom of the device are not blocked.

CAUTION: Take precautions to protect furniture from water damage.

WARNING: Do not use this device without the water chamber installed.

Supplying AC Power to the Device

CAUTION: If this device has been exposed to either very hot or very cold temperatures, allow it to adjust to room temperature (operating temperature) before beginning the following setup procedures. Do not operate the device outside of the operating temperature range shown in the Specifications.

WARNING: Route the wires to avoid tripping.

WARNING: This device is activated when the power cord is connected. Pressing the START/STOP button turns the airflow and humidifier (if activated) on or off.

Complete the following steps to operate the device using AC power.

- 1. Plug the pronged end of the AC power cord into an electrical outlet that is not controlled by a wall switch.
- 2. Plug the power cord's connector into the power inlet on the back of the device.
- 3. Ensure that all connections are secure.

IMPORTANT: To remove AC power, disconnect the power cord from the electrical outlet.

WARNING: Inspect the power cord often for any signs of damage. Discontinue use and replace if damaged.


WARNING: Do not use extension cords with this device.

Setup the Integrated Humidifier

WARNING: Do not use the SleepEasy device without the water chamber in place.

FIRST USE

- 1. Slide the water chamber out from the side of the unit.
- 2. Gently remove the base of the chamber with your hands, being careful not to damage the rubber seal.
- 3. Remove the flapper valve from the air inlet port; unclip the protective cover from the top of the water chamber, as shown in the following diagram. Wash these two parts by hand only in a solution of warm water and mild liquid dishwashing soap. Rinse the parts with clean water and allow them to air dry.



- 4. Wash the water chamber and removable base plate in the dishwasher (top shelf only) or by hand in a solution of warm water and mild liquid dishwashing soap.
- 5. Reassemble. Insert the flapper valve into the air inlet port, and clip the protective cover on top of the water chamber. **Note:** The flapper valve is designed to fit into the chamber in only one orientation. Make sure to place the flapper valve into the chamber in its correct orientation as shown in the above diagram from step 3.
- 6. Fill the chamber to the fill line with distilled water using the funnel as shown here.



CAUTION: Use only distilled water in the chamber. The use of additives or chemicals in the chamber is prohibited.

CAUTION: Do not fill water chamber above the fill line indicated on the side of the water chamber. Damage to the therapy device may occur.

CAUTION: Do not allow the water chamber to sit for any length of time after it has been filled with water. Immediately install the chamber in the device. Allowing water to sit in the chamber (when the chamber is not installed in the device) may cause the chamber to separate from the bottom plate and may cause water leakage.

7. Slide the chamber into place on the side of the unit. Continue to Step 5 of Daily Use.

WARNING: Never touch the heater plate unless the unit is unplugged and the plate has cooled down.

CAUTION: Avoid moving the device when the water chamber has water in it.

DAILY USE

- 1. Slide the water chamber out from the side of the unit.
- 2. Before each use, check the flapper valve operation. Visually inspect the soft flap to make sure it is free to open and close when installed in the water chamber. Check to make sure the soft flap is not pushed back into the flapper valve. Replace the valve if the soft flap is damaged or missing. Make sure the protective cover is clipped into place on top of the water chamber.

Before each use, rinse the chamber with water. Fill the chamber to the fill line with distilled water using the funnel. CAUTION: Use only distilled water in the chamber. The use of additives or chemicals in the chamber is prohibited.

CAUTION: Do not fill water chamber above the fill line indicated on the side of the water chamber. Damage to the therapy device may occur.

CAUTION: Do not allow the water chamber to sit for any length of time after it has been filled with water. Immediately install the chamber in the device. Allowing water to sit in the chamber (when the chamber is not installed in the device) may cause the chamber to separate from the bottom plate and may cause water leakage.

4. Slide the chamber into place on the side of the unit.

WARNING: Never touch the heater plate unless the unit is unplugged and the plate has cooled down. **CAUTION:** Avoid moving the device when the water chamber has water in it.

- 5. Connect the flexible tubing (included with your therapy device system) to the outlet port on the water chamber. **IMPORTANT:** Before each use, examine the flexible tubing for any kinks, damage, or debris. If necessary, clean the tubing to remove the debris. Replace any damaged tubing.
- 6. The ideal humidity setting depends on room temperature and humidity. Initially, a setting of 2 is recommended. You can adjust this setting at any time. Press and hold the humidifier button on the therapy device. The humidifier symbol and setting will appear. Press the + or buttons to change the setting.

IMPORTANT: When the airflow is turned off, the humidifier will automatically shut off. If you restart the airflow, the heated humidification will return to the previous setting. Loss of power to the system will require the user to press the heated humidifier button to reactivate the humidifier.

Display

The display screen is shown here.



The information shown on the display screen is defined as follows:

Icon	Description	Ico	N	Description
35	Indicates that the Blower Hours Time Meter is being displayed.		1	Indicates that the Ramp Start Pressure is being displayed.
X	Indicates that the Therapy Hours Time Meter is being displayed.	cm H ₂		Indicates that a pressure value is being displayed.
>4	When displayed with the Therapy Hours icon, it indicates that the Session Counter is being displayed.	×	•	Indicates that the user may erase the displayed data.
!!!	Indicates that the humidifier is providing heat.	<u>^</u>	7	Indicates that the unit requires user attention.
o	Indicates that the unit is in the Provider Settings Menu.	ے)	Indicates that the Compliance Check Value for the unit is being displayed.
	In Active State, indicates that the Ramp function is in progress. In menus, it is used to indicate Ramp Time Setting alone, or in combination with Ramp Start Pressure.	Ţ	ァ	Indicates that the C-Flex setting is being displayed.

Starting the Device

Note: The numbers shown in the screens throughout this manual are examples only. Actual numbers will vary.

1. Plug the device into an AC power source. The START/STOP, HEAT, RAMP and FLEX buttons light up and the Software Version screen momentarily appears, shown here.



Software Version Screen

2. The next screen to appear is the Calibration screen, shown below. "CAL" will flash on the screen while the device is calibrating its sensors.


Important: Do not connect the breathing circuit until after the unit has finished calibrating.



Calibration Screen

Note: If there was a power failure during a sleep session for longer than 6 hours, the device will recalibrate when the power is restored and the breathing circuit must be removed. If the power failure lasted less than 6 hours during the sleep session, the blower will turn back on at the current settings when power is restored.

3. After successful calibration, the Standby screen will automatically appear, shown below.



Standby Screen

4. Press the START/STOP button () to turn on the airflow. Put on your mask assembly when the air starts to flow.

Note: There will be a short pause after pressing the START/STOP button until the air starts to flow.

5. The Operate screen will then appear, shown here.



Operate Screen

The Operate screen shows the current CPAP setting.

6. Make sure that no air is leaking from your mask into your eyes. If it is, adjust the mask and headgear until the air leak stops. See the instructions provided with your mask for more information.

Note: A small amount of mask leak is normal and acceptable. Correct large mask leaks or eye irritation from an air leak as soon as possible.

7. If you are using the device in bed, try placing the tubing from the device over your headboard. This may reduce tension on the mask.

Note: If you are having trouble with your mask, refer to the instructions supplied with the mask.

Note: You must remove the mask and patient circuit before you get out of bed.

Heated Humidifier Feature

The Integrated Heated Humidifier is a feature that may reduce nasal dryness and irritation by adding moisture and heat to the airflow. This feature can be enabled or disabled by your home care provider.

If the heated humidifier is enabled on your device, press the HEAT button on the top of the device. The Humidifier Setting screen will appear flashing, shown below. You can enter this screen from either the standby or operate screen.

Note: If the heated humidifier feature is disabled by the home care provider, nothing will happen when you press the HEAT button.



Humidifier Setting Screen

You can increase or decrease the humidifier heat settings by pressing the + or - buttons. The possible settings are from 0 (off) to 5 which is the highest setting. Once you have chosen your desired setting, you can either press the HEAT button again or just wait a few moments and the new setting will be saved. The screen will then automatically return to either the standby or operate screen.

Note: This screen will not display if your provider has not enabled the humidifier on your device.

Ramp Feature


You can press the RAMP _____ button to activate the Ramp feature, if your provider has enabled this feature.

The device is equipped with an optional ramp feature that your home care provider can enable or disable. This feature reduces the air pressure when you are trying to fall asleep and then gradually increases (ramps) the pressure until your prescription setting is reached, allowing you to fall asleep more comfortably.

If ramp is enabled on your device, after you turn on the airflow, press the RAMP button on the top of the device. You can use the RAMP button as often as you wish during the night.

Note: If the ramp feature is disabled, nothing will happen when you press the RAMP button.

In order to adjust the Ramp Starting Pressure, press and hold the RAMP \square button for approximately 3 seconds. The Ramp Starting Pressure Screen will appear flashing, shown here. It can be entered from either the standby or operate screen.



Ramp Starting Pressure Screen

You can increase or decrease the ramp starting pressure in 0.5 cm H_2O increments by pressing the + or – buttons. The default setting is 4 cm H_2O . You can adjust the setting from 4 cm H_2O to the CPAP pressure setting. Once you have chosen your desired setting, you can either press the RAMP button again or just wait a few moments and the new setting will be saved. The screen will then automatically return to either the standby or operate screen.

Note: This screen will not display if your provider has not enabled Ramp on your device.

C-Flex Feature

The C-Flex comfort feature allows you to adjust the level of air pressure relief that you feel when you exhale during therapy.

To change this setting, press the FLEX button on the top of the device. The C-Flex Setting screen will appear flashing, shown below. You can enter this screen from either the standby or operate screen.

Note: If the C-Flex feature is disabled by the home care provider, nothing will happen when you press the FLEX button.



C-Flex Setting Screen

You can increase or decrease the setting between 1, 2, or 3 by pressing the + or – buttons. The setting of "1" provides a small amount of pressure relief, with higher numbers providing additional relief. Once you have chosen your desired setting, you can either press the FLEX button again or just wait a few moments and the new setting will be saved. The screen will then automatically return to either the standby or operate screen.

Note: This screen will not display if your provider has not enabled C-Flex on your device.

Navigating the Display Screens

Use the + and - buttons to navigate the display screens.

Note: You can only enter these screens from the standby screen.

Some screens have time-out periods. The screen's timer starts when the screen is initially displayed and is restarted whenever a button is pressed. Unless otherwise specified, all screens timeout after one minute and will return to the standby screen.

Viewing Data on the Patient Data Screens


The following describes the Patient Data screens.



Therapy Usage Hours (Standby Screen)

This screen displays the amount of time that the device provided therapy (with the blower on and the patient connected). The decimal digit is displayed if user hours are less than 10000 so it can be displayed with 0.1 hour resolution. Otherwise, the values between 10000 and 99999 hours can be displayed.


Note: This screen is only for reference. Your home care provider may periodically ask you for this information.



Blower Hours Screen


This screen displays the amount of time that the blower has been active over the life of the device. The decimal digit is displayed if user hours are less than 10000 so it can be displayed with 0.1 hour resolution. Otherwise, the values between 10000 and 99999 hours can be displayed.

Note: This screen is only for reference. Your home care provider may periodically ask you for this information.



Session Counter Screen


This screen displays the number of device therapy sessions that exceeded 4 hours. **Note:** This screen is only for reference. Your home care provider may periodically ask you for this information.



Compliance Check Screen


This screen shows you the 5 digit Compliance Check Value.

Note: Your home care provider may periodically ask you for this information.



Enhanced Compliance VIC: Compliance State Screen

The Enhanced Compliance VIC (Visual Inspection Check) function is used to determine if the device was used for at least 70% of the time over a 30 day period. The Compliance State screen displays a "1" if the device has fulfilled this requirement. **Note:** If a "0" is displayed, the device has not yet met the requirements and the following 4 Enhanced Compliance VIC screens will not be displayed.



Enhanced Compliance VIC: Compliance Month and Day Screen

This screen displays the month and day that the device has met the Enhanced Compliance VIC rules.


Note: Your home care provider may periodically ask you for this information. This screen will not display if the device has not yet met the requirements for Enhanced Compliance VIC.



Enhanced Compliance VIC: Compliance Year Screen

This screen displays the year that the device has met the Enhanced Compliance VIC rules.

Note: Your home care provider may periodically ask you for this information. This screen will not display if the device has not yet met the requirements for Enhanced Compliance VIC.



Enhanced Compliance VIC: Compliance Days > 4 Hours Screen

This screen displays the number of days that the device was used for longer than 4 hours since the device met the 30 days required for the Enhanced Compliance VIC.

Note: Your home care provider may periodically ask you for this information. This screen will not display if the device has not yet met the requirements for Enhanced Compliance VIC.



Enhanced Compliance VIC: Compliance Check Code Screen

This screen shows you the 2 digit Compliance Check Code.

Note: Your home care provider may periodically ask you for this information. This screen will not display if the device has not yet met the requirements for Enhanced Compliance VIC.

System Error Screen

When the unit detects a system error, the System Error screen is displayed as shown below. The blower is turned off and pushbutton functions are disabled. Refer to Alerts and Troubleshooting for more information.



System Error Screen

Device Alerts

High Priority Alert – This alert requires immediate operator response. The alert signal is the backlights on the buttons providing a high priority flashing pattern consisting of a continuous, bright-to-off, two-flash pattern (indicated in the following table as: $\Diamond\Diamond$ $\Diamond\Diamond$ $\Diamond\Diamond$.

ALERT	Visual Indicator	Device Action	Possible Cause	Patient Action
System Error	Backlights: \$\delta \lambda \delta \delta \delta \delta\$ The following symbol displays to indicate that service is required:	The device enters the "Safe state" in which the device power remains on, but the airflow is disabled and the humidifier is turned off.	Device failure	Remove the power supply cord from the device to remove power. Plug the cord back into the device's power inlet to restore power. If the alert continues to occur, contact your home care provider.

Troubleshooting

The table below lists some of the problems you may experience with your device or mask and possible solutions to those problems.

Problem	Why It Happened	What to Do
Nothing happens when you apply power to the device. The backlights on the buttons do not light.	There's no power at the outlet or the device is unplugged.	Check the outlet power and verify that the device is properly plugged in. Make sure the AC power cord is securely connected to the device's power inlet. If the problem continues to occur, contact your home care provider. Return the device to your provider, so they can determine if the problem is with the device.
The device does not operate when you press the () button. The airflow does not turn on.	There may be a problem with the blower.	Make sure the device is powered correctly. If the button backlights turn on when you apply power, but the airflow does not turn on, there may be a problem with your device. Contact your home care provider for assistance. Note: When the device is functioning correctly, after you press the U button, the airflow turns on after a slight delay. This brief delay is
The device's display is erratic.	The device has been dropped or mishandled, or the device is in an area with high Electromagnetic Interference (EMI) emissions.	normal. Unplug the device. Reapply power to the device. If the problem continues, relocate the device to an area with lower EMI emissions (e.g., away from electronic equipment such as cellular phones, cordless phones, computers, TVs, electronic games, hair dryers, etc.). If the problem still occurs, contact your home care provider for assistance.
The Ramp feature does not work when you press the Ramp button.	Your home care provider did not prescribe Ramp for you, or your CPAP pressure is already set to the minimum setting.	If Ramp has not been prescribed for you, discuss this feature with your home care provider to see if they will change your prescription. If your provider has enabled Ramp, but the feature still does not work, check the CPAP setting on your Active Display screen. If CPAP is set to the minimum setting (4.0 cm $\rm H_2O$), or the starting pressure is the same as the prescribed pressure, the Ramp feature will not work.
Air leaks from the vent hole on top of the water chamber.	Improperly assembled flapper (soft) valve. Damaged soft valve.	Remove water chamber and ensure that the flap of the valve operates freely. Reassemble system. Replace soft valve.
You have throat or nose dryness.	The air is too dry.	Increase the humidity setting. When using the heated humidifier feature, refer to these instructions and make sure the humidifier is working properly.
		If the problem still occurs, contact your home care provider for assistance.

Accessories

Contact your home care provider for additional information on the accessories available for your SleepEasy system. When using optional accessories, always follow the instructions enclosed with the accessories.

Adding Supplemental Oxygen

Oxygen may be added at the mask connection. Please note the warnings listed below when using oxygen with the device.

WARNINGS:

- · When using oxygen with this system, the oxygen supply must comply with local regulations for medical oxygen.
- When using oxygen with this system, a Respironics Pressure Valve must be placed in-line with the patient circuit after the oxygen source. Failure to use the pressure valve could result in a fire hazard.
- Oxygen accelerates fires. Keep the device and the oxygen container away from heat, open flames, any oily substance, or other sources of ignition. Do not smoke in the area near the device or the oxygen.
- When using oxygen with this system, turn the device on before turning on the oxygen. Turn the oxygen off before turning the device off. This will prevent oxygen accumulation in the device.
- If administering fixed-flow supplemental oxygen, the oxygen concentration may not be constant. The inspired
 oxygen concentration will vary, depending on the CPAP setting, patient breathing pattern, and leak rate. Substantial
 leaks around the mask may reduce the inspired oxygen concentration to less than the expected concentrations.
 Appropriate patient monitoring should be implemented.
- Do not connect the device to an unregulated or high pressure oxygen source.

Traveling with the System

For your convenience at security stations, there is a note on the bottom of the device stating that it is medical equipment. It may be helpful to bring this manual along with you to help security personnel understand the SleepEasy device.

Note: Before traveling, you must completely empty and dry the water chamber.

Note: The water chamber must be detached from the unit while traveling.

Note: Make sure the soft valve is safely secured to protect from any damage that may occur while traveling.

If you are traveling to a country with a line voltage different than the one you are currently using, a different power cord or an international plug adaptor may be required to make your power cord compatible with the power outlets of the country to which you are traveling. Contact your home care provider for additional information.

Cleaning the Device

WARNING: To avoid electric shock, unplug the device before cleaning it. DO NOT immerse the device in any fluids.

CAUTION: Do not immerse the device in liquid or allow any liquid to enter the enclosure, inlet filter, or any opening.

- 1. Unplug the device, and wipe the outside of the device with a cloth slightly dampened with water and a mild detergent. Let the device dry completely before plugging in the power cord.
- 2. Inspect the device and all circuit parts for damage after cleaning. Replace any damaged parts.

Cleaning or Replacing the Filters

CAUTION: Operating the device with a dirty filter may keep the system from working properly and may damage the device.

Under normal usage, you should clean the gray foam filter at least once every two weeks and replace it with a new one every six months. The white ultra-fine filter is disposable and should be replaced after 30 nights of use or sooner if it appears dirty. DO NOT clean the ultra-fine filter.

CAUTION: Dirty inlet filters may cause high operating temperatures that may affect device performance. Regularly examine the inlet filters as needed for integrity and cleanliness.

- 1. If the device is operating, stop the airflow by pressing the (b) button. Disconnect the device from the power source.
- 2. Remove the filter(s) from the enclosure by gently squeezing the filter in the center and pulling it away from the device.
- 3. Examine the filter(s) for cleanliness and integrity.
- 4. Wash the gray foam filter in warm water with a mild detergent. Rinse thoroughly to remove all detergent residue. Allow the filter to air dry completely before reinstalling it. If the foam filter is torn, replace it. (Only Respironics-supplied filters should be used as replacement filters.)
- 5. If the white ultra-fine filter is dirty or torn, replace it.
- 6. Reinstall the filters, inserting the white ultra-fine filter first if applicable, refer to "Installing the Air Filters".

CAUTION: Never install a wet filter into the device. You must ensure sufficient drying time for the cleaned filter.

Cleaning the Tubing

Clean the tubing daily. Disconnect the flexible tubing from the device. Gently wash the tubing in a solution of warm water and a mild detergent. Rinse thoroughly. Air dry.

Cleaning the Water Chamber

Note: The Flapper (Soft) Valve has a one year lifespan. You must replace this valve after one year.

Note: Hand washing can be performed daily. Dishwashing (water chamber and base plate only) can be performed once a week. The flapper valve and protective cover must be hand washed.

WARNING: Empty and clean the water chamber daily to prevent mold and bacteria growth.

WARNING: Allow the water in the chamber to cool to room temperature before removing the chamber from the humidifier.

- 1. Turn the therapy device off and allow approximately 15 minutes for the heater plate and water to cool.
- 2. Disconnect the tubing from the water chamber.
- 3. Slide the water chamber out of the humidifier platform. Empty any remaining water.
- 4. Gently remove the base of the chamber with your hands, being careful not to damage the rubber seal.
- 5. Remove the flapper valve from the air inlet port; unclip the protective cover from the top of the water chamber. Wash these two parts <u>by hand only</u> in a solution of warm water and mild dishwashing soap. Rinse with clean water and allow them to air dry.
- Wash the water chamber and base plate in the dishwasher (top shelf only) or by hand in a solution of warm water and mild dishwashing soap.
- 7. Inspect all parts for damage prior to reassembly.
- 8. Reassemble the water chamber. Make sure the base plate is fully seated on the water chamber. Also, make sure you install the flapper valve properly into the air inlet port of the water chamber.
 - **Note:** Be sure to reinstall the protective cover on top of the water chamber to prevent debris from entering the water chamber.
- Fill the water chamber to the fill line. Inspect the water chamber for any leaks or damage. Replace the water chamber if any damage is present.

Service

The SleepEasy device does not require routine servicing.

WARNING: If you notice unexplained changes in the performance of this device, if it is making unusual or harsh sounds, if the enclosure is broken, or if water has entered the device, discontinue use, and contact your home care provider.

Specifications

Environmental

	Operating	Storage
Temperature	41 to 95° F (5 to 35° C)	-4 to 140° F (-20 to 60° C)
RELATIVE HUMIDITY	15 to 95% (non-condensing)	15 to 95% (non-condensing)
Atmospheric Pressure	101 to 77 kPa (0 - 2286 m / 0 - 7500 ft)	n/a

Physical

Dimensions: 10° L x 7.5° W x 5.5° H (25 x 19 x 14 cm)

Weight: Approximately 3.7 lbs (1.67 kg)
Water Capacity: 350 ml at recommended water level

Standards Compliance

This device is designed to conform to the following standards:

• IEC 60601-1 General Requirements for Safety of Medical Electrical Equipment

• EN ISO 17510-1 Sleep Apnea Breathing Therapy Devices

• EN ISO 8185 General Requirements for Humidification Systems

Electrical

AC Power Consumption: 100 – 240 VAC, 50/60 Hz, 2.7 A max. Type of Protection Against Electric Shock: Class II Equipment

Degree of Protection Against Electric Shock: Type BF Applied Part

Degree of Protection against Ingress of Water: Device: Drip Proof, IPX1

Mode of Operation: Continuous

Electromagnetic Compatibility: The device meets the requirements of EN 60601-1-2, 2nd edition.

Fuses: There are no user-replaceable fuses.

Humidity

Humidity_{min} Output: 10 mg H₂O/L

Measured @ max flow, 35° C, 15% RH.

Pressure

Pressure Increments: 4.0 to 20.0 cm H₂O (in 0.5 cm H₂O increments)

Pressure Stability: $< 10 \text{ cm H}_2\text{O} (\le 0.5 \text{ cm H}_2\text{O} \text{ peak-to-peak})$

 \geq 10.0 to 20 cm H₂O (\leq 1.0 cm H₂O peak-to-peak)

Measured in accordance with EN ISO 17510-1, 2002 @ 1/3, 2/3, and Pmax with BPM set to 10, 15, and 20 BPM.

Maximum Flow: 35 LPM

Measured in accordance with EN ISO 17510-1, 2002 @ 1/3, 2/3, and Pmax with BPM set to 10, 15, and 20 BPM.

Control Accuracy

 PARAMETER
 RANGE
 ACCURACY

 CPAP
 4 to 20 cm H₂O
 ± 0.5 cm H₂O

Measured at the patient end of the circuit with a Whisper Swivel II exhalation device and no patient flow.

Noise

Sound Pressure Level: < 30 dB(A)

Measured in accordance with EN ISO 17510-1 @ 10 cm H_2O at the patient circuit.

Disposal

Dispose of the device in accordance with local regulations.

How to Contact Respironics

To have your device serviced, contact your home care provider. If you need to contact Respironics directly, call the Respironics Customer Service department at 1-800-345-6443 or 1-724-387-4000.

You can also use the following address:

Respironics 1001 Murry Ridge Lane Murrysville, PA 15668

EMC Information

Guidance and Manufacturer's Declaration - Electromagnetic Emissions - This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Emissions Test	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/Flicker emissions IEC 61000-3-3	Complies	

Guidance and Manufacturer's Declaration - Electromagnetic Immunity – This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

IMMUNITY TEST	IEC 60601 Test Level	Compliance Level	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE	
Electrostatic Discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.	
Electrical fast Transient/burst	±2 kV for power supply lines ±1 kV for input-output lines	±2 kV for supply mains ±1 kV for input/output lines	Mains power quality should be that of a typical home or hospital environment.	
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV for common mode	Mains power quality should be that of a typical home or hospital environment.	
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U _T (>95% dip in U _T) for 0.5 cycle 40% U _T (60% dip in U _T) for 5 cycles 70% U _T (30% dip in U _T) for 25 cycles <5% U _T (>95% dip in U _T) for 5 sec	<5% U _T (>95% dip in U _T) for 0.5 cycle 40% U _T (60% dip in U _T) for 5 cycles 70% U _T (30% dip in U _T) for 25 cycles <5% U _T (>95% dip in U _T) for 5 sec	Mains power quality should be that of a typical home or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.	
Power frequency (50/60 Hz) magnetic field	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical hospital or home environment.	
NOTE: U _T is the a.c. mains voltage prior to application of the test level.				

Guidance and Manufacturer's Declaration - Electromagnetic Immunity - This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Імминіту Теѕт	IEC 60601 Test Level	Compliance Level	ELECTROMAGNETIC ENVIRONMENT -GUIDANCE
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance d = 1.2 √P 80 MHz to 800 MHz d = 2.3 √P 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey*, should be less than the compliance level in each frequency range. ⁵
			Interference may occur in the vicinity of equipment marked with the following symbol: ((**))

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

- a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.
- b Over the frequency range 150 kHz to 80 MHz, the field strengths should be less than 3 V/m.

Recommended Separation Distances between Portable and Mobile RF Communications Equipment and This Device: The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this device as recommended below, according to the maximum output power of the communications equipment.

RATED MAXIMUM POWER OUTPUT OF TRANSMITTER	Separation Distance According to Frequency of Transmitter m			
I KANSMITTEK W	150 кНz то 80 МНz d = 1.2 √Р	80 MHz то 800 MHz d = 1.2 √P	800 МНz то 2.5 GHz d = 2.3 √Р	
0.01	0.12	0.12	0.23	
0.1	0.38	0.38	0.73	
1	1.2	1.2	2.3	
10	3.8	3.8	7.3	
100	12	12	23	

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter; where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1:At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2:These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

Limited Warranty

Respironics warrants that the system shall be free from defects of workmanship and materials and will perform in accordance with the product specifications for a period of two (2) years from the date of sale by Respironics to the dealer. If the product fails to perform in accordance with the product specifications, Respironics will repair or replace – at its option – the defective material or part. Respironics will pay customary freight charges from Respironics to the dealer location only. This warranty does not cover damage caused by accident, misuse, abuse, alteration, and other defects not related to material or workmanship.

Respironics disclaims all liability for economic loss, loss of profits, overhead, or consequential damages which may be claimed to arise from any sale or use of this product. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

This warranty is given in lieu of all other express warranties. In addition, any implied warranties – including any warranty of merchantability or fitness for the particular purpose – are limited to two years. Some states do not allow limitations on how long an implied warranty lasts, so the above limitation may not apply to you. This warranty gives you specific legal rights, and you may also have other rights which vary from state to state.

To exercise your rights under this warranty, contact your local authorized Respironics dealer or contact Respironics at:

1001 Murry Ridge Lane
Murrysville, Pennsylvania 15668-8550
1-724-387-4000

Manufactured for:

Respironics Inc. 1001 Murry Ridge Lane Murrysville, PA 15668 USA

