**PURITAN
BENNETT**

tyco
Healthcare

Sandman Duo™

Sandman DuoST™

Patient Manual

Manual del paciente

Manual do utilizador

Manuel patient

Sandman Duo™, Sandman DuoST™

Patient Manual.....	1
Manual del paciente.....	45
Manual do utilizador.....	89
Manuel patient.....	133

PURITAN
BENNETT

tyco
Healthcare

Sandman Duo™

Sandman DuoST™

Patient Manual

List of Revisions to the Patient Manual Sandman Duo™, Sandman DuoST™

The following pages are included in patient manual
Ref. M-148DFU04-00 revision A.

Pages in process	Revision
Cover	A
List of revisions	A
Table of contents	A
Introduction	A
1 – 44	A

Reference or revision	Description	Date
M-148DFU04-00 revision A.	Manual created	September 2008

Sandman Duo™ and **Sandman DuoST™** are brands owned by Puritan Bennett Inc.

© 2008 Puritan Bennett Inc.

All rights reserved.

TABLE OF CONTENTS

BEFORE YOU START	1
INTRODUCTION	1
SAFETY PRECAUTIONS	1
INDICATIONS FOR USE	2
CONDITIONS OF USE	3
OPERATING PRECAUTION	3
ADVERSE EFFECTS	3
CHECKING THE COMPONENTS AND DEVICE PREPARATION	4
TRANSPORTING THE DEVICE	5
DEVICE DESCRIPTION	6
INSTALLATION	8
DISASSEMBLING THE WATER CHAMBER	8
FILLING THE WATER CHAMBER	9
STANDARD INSTALLATION OF THE DEVICE	10
INSTALLING AN OXYGEN ADAPTER (OPTIONAL)	11
INSTALLATION FOR USE WITH BATTERY POWER OR (OPTIONAL) CIGARETTE-LIGHTER CABLE.....	12
OPERATION.....	13
STARTING TREATMENT	13
STOPPING TREATMENT	15
USING THE RAMP FEATURE	15
DOWNLOADING DATA VIA MEMORY CARD	16
ACCESSING THE DEVICE MENU	17
SANDMAN DUO™ SETTINGS	19
SANDMAN DUOST™ SETTINGS	24
GENERAL INFORMATION	29
REMINDERS ABOUT TO EXPIRE	29
LIST OF AUTHORIZED ACCESSORIES	30
MASKS	30
PATIENT TUBING	30
OTHER ACCESSORIES	30
CLEANING AND MAINTENANCE	32
SANDMAN DUO™ OR SANDMAN DUOST™ DEVICES	32
AIR INTAKE FILTERS	32
CLEANING THE WATER CHAMBER	33
PATIENT TUBING	33
MASK	33
OXYGEN ADAPTER	33
TROUBLESHOOTING	34
HELPFUL HINTS	34
DEVICE MESSAGES	36
TECHNICAL CHARACTERISTICS	37
SANDMAN DUO™ OR SANDMAN DUOST™ DEVICE PERFORMANCE	37
INTEGRATED HUMIDIFIER PERFORMANCE ON SANDMAN DUO™ HC AND SANDMAN DUOST™ HC	37
PHYSICAL CHARACTERISTICS	37
TRANSPORTATION AND STORAGE CONDITIONS	38
OPERATING CONDITIONS FOR THE DEVICE AND HUMIDIFIER	38
DEVICE ELECTRICAL CHARACTERISTICS	38
MEANING OF SYMBOLS DISPLAYED ON THE SCREEN	39
SYMBOLS ON THE DEVICE	40
STANDARDS COMPLIANCE	40

Page intentionally left blank

BEFORE YOU START

Introduction

The **Sandman Duo™** and **Sandman DuoST™** devices are designed for treating patients with Sleep Disordered Breathing including Obstructive Sleep Apnea Syndrome (OSAS). The main characteristic of this pathology is the night-time occurrence of multiple respiratory events, which are brought on by a restriction or obstruction of the upper airways resulting in a disturbance of sleep.

These devices can be operated in either Continuous Positive Airway Pressure (CPAP) mode or Bilevel mode (which delivers two levels of pressure, one inspiratory and the other expiratory) to prevent the upper airways from becoming obstructed and thereby eliminate respiratory arrest. They can be used at home or in a sleep center. The ST model allows a backup frequency to be set.

CAUTION

Federal (U.S.A.) law restricts the Sandman Duo™ and Sandman DuoST™ devices to sale by or on the order of a physician.

Your doctor has decided upon the treatment of your disease. He/she has prescribed the various parameters that the hospital or your home care provider will set on your device. Under no circumstances should you attempt to adjust these prescribed settings.

Please read this manual carefully before using the Sandman Duo™ or Sandman DuoST™ device so that you fully understand the limitations on use of this device.

The information contained in this manual is designed to help you use your device under optimal safety and operating conditions. Contact your medical team or home care provider if you have questions about this equipment or its use.

The device includes an integrated power supply and may be equipped with a heated humidification feature.

Unless otherwise specified, all information contained in this manual applies to Sandman Duo™ and Sandman DuoST™ devices, as well as to Sandman Duo™ HC and Sandman DuoST™ HC, the HC models corresponding to the devices equipped with an integrated humidifier.

Safety precautions

As you read this manual, pay special attention to the special information preceded by the following terms:

WARNING

This signals a risk of death or accident to you or to others.

CAUTION

This indicates that there is a possibility of material damage to this device or any other device.

Note

This indicates a tip for more efficient and practical operation.

WARNINGS

1. Before each use, inspect the Sandman Duo™ or Sandman DuoST™ device following the instructions provided under "Checking the components and device preparation", on page 4.
2. Use only the accessories received from the home care provider. Confirm the 'Directions for Use' instructions for each accessory are in the package and read them carefully.
3. Follow the device installation procedure defined under "Installation", on page 8 of this manual.
4. Keep the Sandman Duo™ or Sandman DuoST™ device far from any source of water. If the device is equipped with an integrated humidifier, fill the water chamber well away from the device to avoid any water spilling onto the device.
5. If supplemental oxygen is to be introduced, carefully follow the instructions and safety precautions for using oxygen, as indicated under "Installing an oxygen adapter (optional)", on page 11.
6. Do not use the Sandman Duo™ or Sandman DuoST™ device near flammable vapors, and in particular do not use the integrated heated humidifier near flammable anesthetics, whether they are present alone or in combination with other gases.
7. Be careful not to obstruct the air outlet or any other opening of the device or respiratory circuit either accidentally or intentionally. Do not introduce liquids or objects into the air outlet.
8. Contact the home care provider if the device is not working properly. Maintenance of this equipment is to be performed by skilled personnel only. Do not try to open the device (risk of electric shock).
9. If the device is equipped with an integrated humidifier, precautions must be taken by the patient when using the water chamber to eliminate the risk of introducing water into the device, which could cause irreversible damage to it. The device must be placed on a horizontal and stable surface, it must not be sloped. The water chamber should always be emptied before moving or transporting the device.
10. If the device is equipped with an integrated humidifier, the humidifier must be cleaned before first use, after any technical intervention, and regularly thereafter, following the instructions provided under "Cleaning and Maintenance", on page 32.
11. Do not use the device if the water chamber is leaking. Leaking may indicate deterioration of the gasket.
12. The water chamber of the integrated humidifier is designed for 'Single Patient Use' only.
13. The heated humidifier must not be used where the ambient temperature is higher than 35 °C (95 °F) because the temperature of the air delivered to the patient could exceed 41 °C (105.8 °F), which could pose a risk of scalding.
14. Do not use this device if the unit appears to be defective or damaged.

Indications for use

The Sandman Duo™ or Sandman DuoST™ device is indicated for the treatment of obstructive sleep apnea in spontaneously breathing patients weighing more than 30kg (66lb) within homecare and hospital environments.

Conditions of use

The Sandman Duo™ or Sandman DuoST™ device is designed for use at home or in sleep centers. It can easily be taken along on trips and used in hotels, and in cars and trucks. When the device is battery-powered, the heated humidification feature cannot be used, but cold humidification is still possible. The water chamber is designed to be filled only with distilled water or tap water that has been boiled and cooled.

Since the Sandman Duo™ and the Sandman DuoST™ are medical electrical devices, please follow the installation instructions contained in this manual concerning electromagnetic compatibility.

Operating precaution

Like all medical electrical devices, the Sandman Duo™ and Sandman DuoST™ are vulnerable to interference from mobile and portable RF communication equipment.

WARNINGS

Like all other devices that generate Continuous Positive Airway Pressure, ensure that the patient is not suffering from:

1. **Pneumoencephallus, trauma or recent surgery with sequela of crano-nasopharyngeal fistula**
2. **Decompensated cardiac insufficiency or hypotension, particularly in case of decreased blood volume or cardiac arrhythmia**

CAUTION

Like all other devices that generate Continuous Positive Airway Pressure, the Sandman Duo or Sandman DuoST device should be used **with caution** in the following conditions:

- Patients who are dehydrated or at risk of dehydration following hydric restriction or diuretic treatment, including voluntary and involuntary changes in treatment.
- Massive epistaxis or history of massive epistaxis
- History of trauma or recent surgery with sequela of crano-nasopharyngeal fistula
- Severe bullous emphysema or emphysema previously complicated by pneumothorax.
- Acute sinusitis, otitis media, or perforated tympanic membrane
- Severe claustrophobia.
- The patient must not use the heated humidification feature if the patient presents a deviation of the supraglottic airways.

Adverse effects

Contact your medical team if the following symptoms appear while using the Sandman Duo™ or Sandman DuoST™ device: airway or nasal passage dryness, skin sensitivity, runny nose, earache, sinus discomfort, daytime drowsiness, mood changes, disorientation, irritability or memory loss and/or epistaxis.

Checking the components and device preparation

The Sandman Duo™ or Sandman DuoST™ device is delivered in a carrying case packaged with a detachable AC power cord, patient tubing, a filter, a patient manual, and a memory card. Check the integrity of packaging when you receive the device.

Carefully examine the Sandman Duo™ or Sandman DuoST™ device and each accessory. Check that the patient circuit is undamaged and that the mask supplied by the home care provider has no cracks or other defect in the shell or bellows chamber.

If the device is equipped with an integrated humidifier, the water chamber is delivered inside the carrying case, separate from the device. Check the condition of the water chamber. If it is damaged, contact your home care provider.

To fill the water chamber and attach to the device, follow the instructions provided in steps 1 to 4, "" (page 44).

If the device is not equipped with an integrated humidifier, the cover is delivered inside the carrying case, separate from the device. Check the condition of the cover.

To install the cover on the device, set the base of cover in the dummy plate and tip it forward to lock it in place.


Figure 1 - Installing the device cover

If supplemental oxygen is to be introduced, you will also need an oxygen adapter (optional) and a special one-way valve to prevent oxygen from accumulating in the device when the device is halted (see "", page 44)".

Inspect the Sandman Duo™ or Sandman DuoST™ device and the detachable AC power cord to make sure they are free of visible defects. Plug one end of the detachable AC power cord into the power port of the device (item 6 of Figure 4 page 7) and the other end into a wall power outlet, as described under "" (page 44). The device will be in the same mode it was in (on/off) when last disconnected from the power cord.

In standby mode, the time and date are displayed (see Figure 2).

Figure 2 – View of display in standby mode

Push the on/off button (Item 4 of Figure 3 on page 7) to start the device. The screen illuminates. On the Sandman Duo™ or Sandman DuoST™ devices, the estimated mask pressure is displayed. The backlight for the display will come on and remain illuminated

for 30 seconds after the last key press. The amount of time for the backlight on the display is not adjustable.

Press and hold the on/off button for several seconds to stop the device and unplug the power cord from the power outlet.

WARNING

If any anomaly is observed in the device (noisy, unexpected display, etc.), do not begin treatment, but contact the home care provider.

Transporting the device

If the device is equipped with an integrated humidifier, you must detach the water chamber, empty it completely, and pack the chamber in the carrying case separate from the device.

CAUTION

If water enters the device it can cause irreversible damage. To eliminate any risk of introducing water into the device, the water chamber needs to be emptied before moving or transporting the device. Be careful not to tilt the device when removing filled water chamber.

To detach the water chamber from the device, follow the instructions provided in Steps 1. to 3., "Disassembling the water chamber" (page 8). Then empty the water chamber completely.

Next, pack the water chamber and the device in the compartments provided for them in the carrying case.

If the device is not equipped with an integrated humidifier, detach the device cover and pack it in the carrying case separate from the device.

To detach the cover from the device, place your fingers on the grip marks on either side of the cover and press down to separate the cover as shown in Figure 1. Remove the cover.

Next, pack the cover and the device in the compartments provided for them in the carrying case.

DEVICE DESCRIPTION

The **Sandman Duo™ or Sandman DuoST™ device** is described in the following. It is equipped with custom accessories and includes an integrated electrical power supply. See Figures 3 to 7 on page 7.

1. Display

This allows information such as the pressure, device settings and the time to be viewed. For further information, refer to "", on page 44.

2. Left control button

This allows the selection of the feature offered at the lower left of the display. For further details, refer to "", on page 44.

3. Ramp button

This allows the ramp time setting to be accessed, but only when the device is in standby mode (for further information, see "", on page 44). It is also used to select parameters from the device-setting menus.

4. On/off button

This allows the device to be turned on or off when it is in standby mode. It is used also to exit the device-setting menus.

5. Right control button

This allows the selection of the feature offered at the lower right of the display. For further details, refer to "", on page 44.

6. Power port

This is where the detachable power cord is connected. For further information, refer to "", on page 44.

7. External battery port

This allows the device to be powered by an external battery or a cigarette-lighter power cable. For further information, refer to "Installation for use with battery power or (optional) cigarette-lighter cable", on page 12.

8. Air intake filters

These filters prevent dust in the airflow from entering the device. Refer to "Cleaning and Maintenance", on page 32, for further information about how to maintain these filters.

9. Heating element (depending on model)

The built-in heated humidifier located beneath the water chamber allows the water in the water chamber to be heated.

WARNING

Avoid all contact with this element.

10. Serial cable/USB connector

This connector available on Sandman Duo™ and Sandman DuoST™ devices is for the use of your physician or home care provider. Do not connect anything to it.

11. Memory card connector

This connector allows the insertion of the memory card to download the data in the device's memory. For further information, refer to "", on page 44.

12. Outlet connector

Located on the cover of the device or above the water chamber, this is the air outlet to which the patient tubing will be connected. For further information, refer to "Installation", on page 8.

13. Water chamber (depending on model)

This is a transparent plastic container located at the rear of the device. It is designed to hold the water needed for humidification and is marked to show the maximum fill level for water.

The components in the following figures are described on page 6.


Figure 3 – Front view


Figure 4 – View of the right side


Figure 5 – View of the rear side


Figure 6 – Rear view of a device without water chamber


Figure 7 – Rear view of a device with water chamber

INSTALLATION

This section indicates how to connect the Sandman Duo™ or Sandman DuoST™ accessories. Before proceeding with installation, inspect the Sandman Duo™ or Sandman DuoST™ device and its accessories as described under "Checking the components and device preparation", on page 4.

WARNINGS

- **Place the Sandman Duo™ or Sandman DuoST™ device on a horizontal stable surface in a clean environment.**
- **Avoid placing the device too close to a wall, which could interfere with airflow into the device.**

Disassembling the water chamber

If the device is not equipped with an integrated humidifier, skip to paragraph 1.b of "", on page 44.

If the Sandman Duo™ or Sandman DuoST™ device has an integrated humidifier (HC model), follow this procedure:

1. First, turn off the device by holding down the on/off button for several seconds and remove the patient circuit if necessary. As the integrated humidifier is located at the rear of the device, place your fingers on the grip marks on either side of the cover and press in to separate the water chamber.


Figure 8 – Separating the water chamber

2. Tip the water chamber back as shown in Figure 9. Remove the water chamber from the machine.


Figure 9 – Tipping the water chamber

3. Unlock the water chamber cover by pressing on the two clips as shown (Figure 10a) and remove the water chamber cover from the base (Figure 10b). The silicone gasket may be left in place on the base. Be careful since the water chamber may be hot.


Figure 10 – Opening the water chamber

Filling the water chamber

CAUTION

Before using the water chamber, check that it is clean. If it is not, refer to "Cleaning and Maintenance", on page 32 of this manual.

1. Carefully pour water into the water chamber up to the maximum fill level.

WARNING

Fill the water chamber away from the device to avoid spilling water on the device.

CAUTION

- **Be careful not to exceed the maximum fill level.**
- **Either distilled water or boiled and cooled tap water may be used in the water chamber.**
- **Do not add any other product to the water in the water chamber.**
- **Do not use alkaline solution (physiological saline solution).**
- **Use only room temperature water; do not fill the water chamber with hot or chilled water.**

2. If you have removed the silicone gasket, replace it securely on the water chamber base within the channel provided for this purpose (Figure 11). Otherwise, the device will not deliver air correctly.


Figure 11 – Replacing the silicone gasket

3. Position the rear edge of the water chamber cover on the water chamber hinge and reseat the upper part of the water chamber on the base. Lock the water chamber.
4. Attach the water chamber base to the heating element and tip it forward to lock it in place, being careful not to spill any water (Figure 12).


Figure 12 – Replacing the water chamber

WARNINGS

- **There is a risk of air leakage if the water chamber or the cover of the Sandman Duo™ or Sandman DuoST™ is not reassembled correctly, which could cause the applied pressure to deviate from the prescribed pressure.**
- **The humidifier has a plate that, like the base of the water chamber, may become hot during normal operation and should not be touched.**

Standard installation of the device

To install the device, follow these steps:

1. Connecting the patient circuit
 - a. If the Sandman Duo™ or Sandman DuoST™ device is equipped with an integrated humidifier, once the water chamber is filled, connect one end of the patient circuit to the outlet connector located on the water chamber of the humidifier, as shown in Figure 13a.
 - b. If the Sandman Duo™ or Sandman DuoST™ is not equipped with an integrated humidifier, connect one end of the patient circuit to the outlet connector located at the rear of the device, as shown in Figure 13b.


Figure 13 – Connecting the patient circuit to the device

WARNING

Do not block the air intake.

2. Prepare the mask according to the instructions provided in the operating instructions that came with the mask.

Note: Sandman Duo™ or Sandman DuoST™ can be used with full face masks recommended by Puritan Bennett.

3. Connect the mask to the other end of the patient circuit.
4. Plug the detachable power cord into the power port of the device (Item 6 of Figure 4 on page 7).

WARNING

Before connecting the power cord to the device, make sure that the power cord is undamaged.

5. Plug the other end of the detachable cord to into the power outlet. The display illuminates showing the time and date as soon as the connection is made.
6. The device is ready to operate. Refer to "" on page 44.

Installing an oxygen adapter (optional)

If supplemental oxygen is required, follow the manufacturer's installation instructions of the special one-way valve to prevent oxygen from accumulating in the device when the device is halted, and then replace Step 3 of the standard installation procedure with the following step:

3. Connect one end of the oxygen adapter to the free end of the patient circuit and the other end to the mask.

Continue with steps 4 through 6 of the standard installation procedure.

WARNINGS

- **When using oxygen, always follow the instructions of the medical team or home care provider.**
- **Do not smoke in the presence of oxygen.**
- **Do not inject oxygen into the device's air intake.**
- **Carefully follow the instructions for starting and stopping treatment, defined under " on page 44.**
- **If you use an oxygen concentrator or liquid oxygen unit, stop the flow of oxygen when the device is not operating. If the oxygen concentrator remains on when the device is turned off, the oxygen delivered into the patient circuit could accumulate inside the device, creating a fire hazard.**

CAUTION

For a given oxygen flow, the concentration of oxygen inhaled varies as a feature of the pressure settings, your breathing, the type of mask used and the leak rate. This precaution applies to most Continuous Positive Airway Pressure devices.

Installation for use with battery power or (optional) cigarette-lighter cable

The Sandman Duo™ or Sandman DuoST™ device can be powered by a 12-volt DC battery using an optional cable specially designed for this purpose.

CAUTION

- **Make sure you use a battery power cable that is compatible with Sandman Duo™ or Sandman DuoST™ device. Using a non-approved cable will increase the risk of damaging both the device and the battery.**
- **Use only a 12-volt battery power source and observe the correct connection polarity.**
- **Connecting a mains voltage DC power supply using 12-volt input is not allowed.**

Note: When the device is battery-powered, the heated humidification feature cannot be used, but a cold humidification is possible.

To install the battery cable, replace Steps 4 and 5 of the standard installation procedure with the following steps:

4. Plug the battery power cable into the external battery port of the device (Item 7 of Figure 4 on page 7).
5. Plug the other end of the cable directly into the battery.

The device can be connected simultaneously to an electrical outlet and a backup battery, in which case follow Steps 4 and 5 of the standard installation procedure to connect to electrical power and the steps described above to connect to battery power.

The device will then operate on electrical power and in the event of a power failure the backup battery will automatically take over without interrupting device operation.

Note: It is recommended that you use a battery that is fully charged. The Sandman Duo™ or Sandman DuoST™ device does not recharge a battery when connected to electrical power. Check the battery charge before each use.

The Sandman Duo™ or Sandman DuoST™ device also can also be powered from a cigarette-lighter socket using an optional cable specially designed for this purpose.

For this configuration, replace Steps 4 and 5 of the standard installation procedure with the following steps:

4. Plug the cigarette-lighter cable into the external battery port of the device (Item 7 of Figure 4 on page 7).
5. Plug the other end of the cable directly into the cigarette-lighter socket.

OPERATION

This section describes how to use the Sandman Duo™ or Sandman device when it is installed as directed in the preceding section. Included are the instructions for how to use the pressure ramp feature, how to access settings information to view it or change it, and how to download data in memory to a memory card.

WARNING

In the "Starting Treatment" and "" sections, it is essential that the steps are followed in the order given. Failure to do so could cause an accident, particularly in cases where oxygen is used.

Starting Treatment

1. To start the device, press and release the on/off button (item 4 of Figure 3 on page 7). The screen illuminates and displays the estimated mask pressure with the Sandman Duo™ or Sandman DuoST™ device (see Figure 14).


Figure 14 – View of the display when the device is operating

Note: The unit of measure displayed is set by your home care provider at hPa or cmH₂O (as above).

The symbols that may be displayed on the screen are summarized in the table entitled "Meaning of symbols displayed on the screen", on page 39.

If the Sandman Duo™ or Sandman DuoST device is equipped with an integrated heated humidifier, the heating feature is automatically enabled when the humidification level setting is set between 1 and 10. The heated humidifier indicator then displays on the screen. When the heating feature is disabled, the humidification symbol displays without steam . To adjust humidification, refer to "Accessing the device menu", on page 17.

Note: On a setting of ten (10) with the chamber filled to the maximum fill line, the water will last up to 8 hours. The length of time the water will last is dependant on several factors such as environment, leak rate, patient's breathing pattern, etc.

WARNINGS

When using a heated humidifier:

- **Check that the water chamber is properly seated. If it is not, re-seat it, because contact with the heating element could cause scalding.**
- **Check that the water chamber is filled to the maximum fill level to ensure humidification throughout the entire period of use.**
- **The water chamber is designed to be 'single patient use' only.**
- **Do not tilt the Sandman Duo™ or Sandman DuoST™ device when using the humidifier.**

CAUTION

- **If symptoms of nasal dryness or dry mouth appear, gradually raise the humidity until the symptoms disappear.**
 - **If the condensation in the water chamber or the patient circuit appears to be excessive, or if water droplets appear in the mask, select a lower humidity setting (see "", on page 44).**
 - **The water chamber needs to be emptied before moving or transporting the device. This eliminates the risk of introducing water into the device, which could cause irreversible damage.**
2. Place the mask on your head according to the instructions for use that came with the mask. If the mask has an exhalation vent, it has an opening through which exhaled air is expelled so that the air cannot be re-breathed.

Breathe normally.

If the message "MASK UNPLUGGED" displays on the screen, this means that the mask is not connected properly. Reposition your mask correctly to eliminate the leaks, and press any of the four buttons on the front of the device. The device will restore the pressure set and the message "MASK UNPLUGGED" will disappear.

WARNINGS

- **Never block the mask exhalation vent, which allows air to be expelled continually. When the device is operating, the air produced flushes the exhaled air out through the mask vent hole. When the device is not operating, insufficient fresh air is produced in the mask and there is a risk of re-breathing exhaled air, which under certain circumstances could give rise to a suffocation event within several minutes.**
- **In the first default case, if the patient is equipped with a nasal mask, the device resistance is low enough that the patient may exhale through the machine or the patient may even open his or her mouth. If the patient is equipped with a full face mask, the mask must have an anti-asphyxia valve.**
- **If upper airway dryness is experienced and the device is not equipped with an integrated humidifier, ask your home care provider to supply you a device with integrated humidifier.**

3. If oxygen is used, start the flow.

WARNING

It is essential that the device be turned on and generates air flow before the oxygen flow is started so that oxygen does not accumulate in the device.

Stopping treatment

When the Sandman Duo™ or Sandman DuoST™ device is on, take the following steps to stop treatment:

1. If oxygen is used, stop the flow.

WARNING

You must stop the oxygen flow before turning off the device so that oxygen does not accumulate in the device.

2. Remove the mask.
3. Press and hold the on/off button for several seconds (Item 4 of Figure 3 on page 7). The device will go into standby mode and the humidification feature will stop automatically.

Using the ramp feature

The Ramp feature allows for a gradual rise in pressure to help you go to sleep.

In CPAP mode, treatment begins at a low pressure, called the comfort pressure, and rises from there to the pressure prescribed by the physician.

In Bilevel mode, the patient pressure is reduced to the comfort pressure and then raised to the inspiratory positive airway pressure (IPAP) pressure and the expiratory positive airway pressure (EPAP).

The ramp feature will start automatically once the device has been operating for 5 minutes as long as the ramp time is not zero. The symbol  will start blinking on the left of the display at startup and will become steady once the ramp is active.

The ramp time may be set from 0 up to the maximum prescribed ramp time when the device is in standby. To do this, follow the setting instructions for the Sandman Duo™ (page 19) and the Sandman DuoST™ (page 24).

The ramp time setting may be accessed directly by pressing the ramp button on the front of the device (Item 3 of Figure 3 on page 7), but only when the device is in standby. If the ramp time is set to 0, the ramp feature is disabled. To reactivate the ramp feature, stop and restart the device by pressing twice the on/off button (Item 4 of Figure 3 on page 7). This feature can be used if your sleep is interrupted during the night and you have difficulty getting back to sleep.

Note: If the ramp time is set at zero, the pressure will be the same as the set pressure when the device is operating in CPAP mode or the same as the inspiratory positive airway pressure (IPAP) or the expiratory positive airway pressure (EPAP) when the device is operating in Bilevel mode.

The maximum ramp time is prescribed by your medical team and can be set by your home care provider in 5-minute increments from 0 (no gradual rise in pressure) to 45 minutes. If 0 was prescribed, the ramp function is disabled.

Similarly, the comfort pressure can be set from 3 cmH₂O to the prescribed pressure when in CPAP mode, and from 3 cmH₂O to the expiratory positive airway pressure (EPAP) when in Bilevel mode. For detailed instructions, see pages relating to settings: page 19 for Sandman Duo™ and page 24 for Sandman DuoST™.

Downloading data via memory card

An external memory card can be used to download the most recent compliance data in the device memory. Please refer to the instructions given by your home care provider to use this card.

Note: The Downloading data feature can be accessed only when the device is in standby mode.

If the device is running, stop it by pressing the on/off button for several seconds.

Insert the memory card supplied by the home care provider in the memory card expansion slot (item 11 of Figure 5 on page 7). The memory card symbol  appears in the upper right corner of the screen. If the symbol is blinking, this means that the card is not working and needs to be returned to the home care provider.

Data download begins automatically. An indicator allows the progress of the data downloaded to be monitored. Once the data has been recorded, a message displays on the device "DOWNLOAD FINISHED. REMOVE MEMORY CARD". The memory card can then be given to the home care provider.

Note: If the card is left inside the device, the machine cannot be turned on and no feature can be accessed.

CAUTION

Use memory cards referenced in the list of authorized accessories, on page 30.

Accessing the device menu

The display allows the **Sandman Duo™** or **Sandman DuoST™** information (Item 1 of Figure 3 on page 7) to be viewed.

The parameters are accessible only when the machine is in standby or operating mode. The parameters that are accessible are grouped on several screen pages indicating:

- device settings related to prescribed treatment
- general device settings, such as date and time
- compliance data
- patient reminders - The device has a patient reminder feature to remind patients of key events such as when to change the filter or mask, how to control the device, and how/when to transfer compliance data to the memory card.

Description of a screen page

As shown in Figure 15, each screen page includes:

- a first line or page header
- up to 4 parameter lines with associated values
- a final line showing the various symbols used for such features as changing a parameter value, moving about a page, or moving from page to page (see table entitled "Meaning of symbols displayed on the screen", on page 39).


Figure 15 – Diagram of screen page and buttons on device front

Description of buttons and symbols

The four buttons on the front of the device allow the following:

- access the device menu
- scroll through various pages
- change the value of certain parameters on those pages.

The closed padlock symbol at the bottom of the screen where the page header is selected indicates that the values displayed on the page cannot be modified. If this symbol is displayed at the bottom of the screen where a parameter is selected, the value of this parameter cannot be modified.

The open padlock symbol at the bottom of the screen where the page header is selected shows that the page displays settings. If this symbol is displayed at the bottom of the screen where a parameter is selected, the value of this parameter can be modified.

Using the **left and right control buttons**, any of the symbols offered opposite the last line of the displayed page is selectable:

- The General Info menu can be accessed by pressing the left control button .
- The settings menu is accessed by pressing the right control button opposite the settings symbol .
- When the page header is selected (i.e., when the background is shaded), the control buttons can be used to change the page. Press the right control button  to access the next page, and press the left control button  to go back to the preceding page.
- When a parameter is selected (i.e., when its background is shaded), its value can be increased by pressing the right control button  or decreased by pressing the left control button . To scroll through the available values for the displayed parameter, hold down either of these buttons.

The next parameter on the page is accessible by pressing the **ramp button** . This button also returns you to the top of the page when the last parameter on the page is selected.

Pressing the **on/off button**  allows you to return to the top of the page, and pressing it a second time allows you to exit the menu.

If no button is touched for 30 seconds, the device exits the menu and the screen displays the pressure or the time.

This section provides a general description of the settings that are common to all Sandman Duo™ and Sandman DuoST™ devices.

For settings specific to the **Sandman Duo™** device, see "Sandman Duo™ settings", on page 19.

For settings specific to the **Sandman DuoST™** device, see "Sandman DuoST™ settings", on page 24.

Refer also to "General Information" section for the Sandman Duo™ device and for the Sandman DuoST™ device on page 29.

The patient reminder feature is intended to help you remember key events, such as filter and mask changes, device inspections, and memory card compliance data downloads. This feature, which is identical on both Sandman Duo™ and Sandman DuoST™ devices, is described under "Reminders about to expire", on page 29.

Sandman Duo™ settings


See “” on page 44 for general descriptions of a screen page and the buttons and symbols.

Note: In the following screens, the values displayed are given by way of example.

Pressing the right control button  (Item 5 of Figure 3 on page 7) causes the following screen to display.

1) Humidification level

Views of the display

Note: This screen displays only if the device is equipped with the humidification feature.

The  symbol indicates that the heat level can be modified. The heat level can be set in increments of 1 from 0 to 10.

Press the ramp button to select the setting.

The value may be increased by using the right control button  and decreased by using the left control button . When the heat setting is at 10, the heating element delivers maximum heat.

Once the settings are made, select the page header by using the ramp button or the on/off button, then press the right control button  to access the next page.

Note: If the humidifier indicator at the top of the screen blinks (without the humidification level value) when the device is in standby (i.e., when the time is displayed) or when the device is operating (i.e., when the pressure is displayed),

pressing on the settings button  will cause the blinking  symbol to appear. If the device is connected to a battery, the heated humidification function cannot be used and the symbol blinks. In the other cases, contact your home care provider.

2) Settings

In CPAP mode

View of the display

SETTINGS	CPAP
COMFORT PRESSURE	6.0 cmH ₂ O
RAMP TIME	15 min
PRESSURE	9.0 cmH ₂ O

Comfort pressure: Level of pressure produced by the device when the ramp feature starts up.

Ramp time: Time the device takes to reach the prescribed pressure when starting from the comfort pressure.

Pressure: Set pressure, level of pressure prescribed, which may be altered by $\pm 1 \text{ cmH}_2\text{O}$. The set pressure will display only if the Pressure Alter feature has been selected.

Note: This page will display only if the maximum ramp time is other than 0 (see Screen Parameters) and/or if the Pressure Alter function is activated. (see paragraph “3) Parameters”).

The first parameter is accessible by pressing the ramp button.

The page will change as follows:

Views of the display

SETTINGS	CPAP
COMFORT PRESSURE	6.0 cmH ₂ O
RAMP TIME	15 min
PRESSURE	9.0 cmH ₂ O
-	+

The comfort pressure value is selected and can be set in increments of 0.5 cmH₂O from 3 cmH₂O up to the set pressure, but only when the device is in standby. Press the right control button **+** to raise the setting and the left control button **-** to lower it.

Press the ramp button to access the next parameter on the menu.

SETTINGS	CPAP
COMFORT PRESSURE	6.0 cmH ₂ O
RAMP TIME	15 min
PRESSURE	9.0 cmH ₂ O
-	+

The ramp time value can be set in increments of 5 minutes from 0 (ramp disabled) to the maximum ramp time prescribed, but only when the device is in standby. Press the right control button **+** to raise the setting and the left control button **-** to lower it.

Note: This setting can also be accessed by pressing the ramp button, but only when the device is in standby.

Press the ramp button to access the next parameter on the menu.

SETTINGS	CPAP
COMFORT PRESSURE	6.0 cmH ₂ O
RAMP TIME	15 min
PRESSURE	9.0 cmH ₂ O
-	+

The set pressure value does not display unless the Pressure Alter feature is selected. It can be set in increments of 0.5 cmH₂O to the prescribed pressure \pm 1 cmH₂O. Press the right control button **+** to raise the value and the left control button **-** to lower it within the limit from 3 cmH₂O to 20 cmH₂O.

In Bi-level mode

Views of the display

SETTINGS	BILEVEL
COMFORT PRESSURE	6.0 cmH ₂ O
RAMP TIME	15 min
---	---

This page displays only when at least one parameter is adjustable. Otherwise, it does not display.

Comfort pressure: Level of pressure produced by the device when the ramp feature starts up.

Ramp time: Set time to reach the inspiratory positive airway pressure (IPAP) and the expiratory positive airway pressure (EPAP) when starting from the comfort pressure.

Note: This page does not display unless the maximum ramp time is other than 0 (see the Parameters screen). The Pressure Alter feature is unavailable in Bilevel mode.

The first parameter on the menu can be accessed by pressing the ramp button. The symbols displayed in the lower part of the screen will change as shown in the following.

SETTINGS	BILEVEL
COMFORT PRESSURE	6.0 cmH ₂ O
RAMP TIME	15 min
-	+

The comfort pressure value displays only if the ramp can be activated. The value selected can be set in increments of 0.5 cmH₂O from 3 cmH₂O up to the expiratory positive airway pressure (EPAP), but only when the device is in standby. Press the right control button **+** to raise the setting and the left control button **-** to lower it.

Press the ramp button to access the next parameter on the menu.

SETTINGS	BILEVEL
COMFORT PRESSURE	6.0 cmH ₂ O
RAMP TIME	15 min
-	+

The ramp time value can be set in increments of 5 minutes from 0 (ramp disabled) to the maximum ramp time prescribed, but only when the device is in standby. Press the right control button **+** to raise the setting and the left control button **-** to lower it.

Note: This setting can also be accessed by pressing the ramp button, but only when the device is in standby.

Once the settings are made, select the page header by using the ramp button or the on/off button, then press the right control button **→** to access the next page.

3) Parameters

Note: The  symbol indicates that the values displayed on the page cannot be modified.

In CPAP mode

View of the display

PARAMETERS	CPAP
PREScribed P.	8.0 cmH ₂ O
PRESSURE ALTER	<input checked="" type="checkbox"/>
MAX. RAMP TIME	45 min

Prescribed P.: Level of pressure prescribed by your physician.

Pressure Alter: Feature enabled by the physician or home care provider in order to change the prescribed pressure by ± 1 cmH₂O. The symbol indicates that the Pressure Alter feature is enabled, and the symbol indicates that it is disabled.

Max ramp time: Maximum ramp time of from 0 to 45 minutes set by the home care provider.

In Bilevel mode

Views of the display

PARAMETERS	BILEVEL
MAX. RAMP TIME	45 min
EPAP	5.0 cmH ₂ O
IPAP	15.0 cmH ₂ O

Max ramp time: Maximum ramp time of from 0 to 45 minutes set by the home care provider.

EPAP: Set expiratory pressure (low-pressure setting for expiratory positive airway pressure).

IPAP: Set inspiratory pressure (high-pressure setting for inspiratory positive airway pressure).

Press the right control button  to access the next page.

Press. Rise Time: Pressure rise time set by physician. Time the device takes to rise from the expiratory positive airway pressure (EPAP) to the inspiratory positive airway pressure (IPAP).

Inspi. Sens.: Inspiratory sensitivity set by physician (detection sensitivity of the start of the inspiratory phase of respiration).

Expi. Sens.: Expiratory sensitivity set by physician (detection sensitivity of the start of the expiratory phase of respiration).

Press the right control button  to access the next page.

4) Patient reminders

This screen displays only when the physician or home care provider has set at least one frequency for one type of reminder. Then, the reminder activation date is shown opposite the corresponding reminder.

The display is the same whether the device is set to CPAP mode or Bilevel mode.

Note: The  symbol indicates that the values displayed on the page cannot be modified.

View of the display

REMINDER	
FILTER	09.07.2008
MASK	15.11.2008
SERVICE	09.06.2009
MEMORY CARD	15.08.2008

Filter: Provides a reminder of the date when the filter should be changed.

Mask: Provides a reminder of the date when the mask should be changed.

Service: Provides a reminder of the date when the device pressure should be inspected.

Memory card: Provides a reminder of the date when the memory card should be returned to the home care provider after retrieving the data.

Press the right control button  to access the next page.

5) Setting the clock

The display is the same whether the device is set to CPAP mode or Bilevel mode.

Views of the display

CLOCK	
HOUR	13
MINUTE	45
← →	□

The time can be set when the device is in standby only.

Press the ramp button to select the hour.

CLOCK		
HOUR	13	
MINUTE	45	
-	□	+

The value is based on 24 hours and is set in increments of 1. Press the right control button **+** to raise the setting and the left control button **-** to lower it.

Press the ramp button to select the minutes.

CLOCK		
HOUR	13	
MINUTE	45	
-	□	+

The minutes value is based on 60 minutes and is set in increments of 1. Press the right control button **+** to raise the number of minutes and the left control button **-** to lower it.

*Once the settings are made, select the page header by using the ramp button or the on/off button, then press the right control button **→** to access the next page.*

6) Setting the date

The display is the same whether the device is set to CPAP mode or Bilevel mode.

Views of the display

CALENDAR	
DATE	18
MONTH	9
YEAR	2008
← →	□

The time can be set when the device is in standby only. Press the ramp button to select the day.

CALENDAR		
DATE	18	
MONTH	9	
YEAR	2008	
-	□	+

Press the right control button **+** to raise the day setting and the left control button **-** to lower it.

Press the ramp button to select the month.

CALENDAR		
DATE	18	
MONTH	9	
YEAR	2008	
-	□	+

Press the right control button **+** to raise the month setting and the left control button **-** to lower it.

Press the ramp button to select the year.

CALENDAR		
DATE	18	
MONTH	9	
YEAR	2008	
-	□	+


Press the right control button **+** to raise the year setting and the left control button **-** to lower it.

*Once the settings are made, select the page header by using the ramp button or the on/off button, then press the right control button **→** to access the next page.*

7) Setting the contrast

The display is the same whether the device is set to CPAP mode or Bilevel mode.

Views of the display

Press the ramp button to select the value.

The contrast is set at 5 by default, but can be set from 0 to 10 in increments of 1. Increase the value by using the right control button **+** or decrease it by using the left control button **-**. When the setting is at 10, the contrast is the maximum.

Once your settings are made, select the page header by using the ramp button or the on/off button, then press the on/off button to exit the menu.

Sandman DuoST™ settings


See "" on page 44 for general descriptions of a screen page and the buttons and symbols.

Note: In the following screens, the values displayed are given by way of example.

Pressing the right control button  (Item 5 of Figure 3 on page 7) causes the following screen to display.

1) Humidification level

Views of the display

Note: This screen displays only if the device is equipped with the humidification feature.

The  symbol indicates that the heat level can be modified. The heat level can be set in increments of 1 from 0 to 10.

Press the ramp button to select the setting.

The value may be increased by using the right control button  and decreased by using the left control button . When the heat setting is at 10, the heating element delivers maximum heat.

Once the settings are made, select the page header by using the ramp button or the on/off button, then press the right control button  to access the next page.

Note: If the humidifier indicator at the top of the screen blinks (without the humidification level value) when the device is in standby (i.e., when the time is displayed) or when the device is operating (i.e., when the pressure is displayed),

pressing on the settings button  will cause the blinking  symbol to appear. Contact your home care provider. If the device is connected to a battery, the heated humidification function cannot be used and the symbol blinks.

2) Settings

In CPAP mode

View of the display

SETTINGS	CPAP
COMFORT PRESSURE	6.0 cmH ₂ O
RAMP TIME	15 min
PRESSURE	9.0 cmH ₂ O

Comfort pressure: Level of pressure produced by the device when the ramp feature starts up.

Ramp time: Time the device takes to reach the prescribed pressure when starting from the comfort pressure.

Pressure: Set pressure, level of pressure prescribed, which may be altered by ± 1 cmH₂O. The set pressure will display only if the Pressure Alter feature has been selected.

Note: This page will display only if the maximum ramp time is other than 0 (see Screen Parameters) and/or if the Pressure Alter function is activated. (see paragraph "3) Parameters").

The first parameter is accessible by pressing the ramp button.

The page settings can be changed as follows:

Views of the display

SETTINGS	CPAP
COMFORT PRESSURE	6.0 cmH ₂ O
RAMP TIME	15 min
PRESSURE	9.0 cmH ₂ O

The comfort pressure value is selected and can be set in increments of 0.5 cmH₂O from 3 cmH₂O up to the set pressure, but only when the device is in standby. Press the right control button to raise the setting and the left control button to lower it.

Press the ramp button to access the next parameter on the menu.

SETTINGS	CPAP
COMFORT PRESSURE	6.0 cmH ₂ O
RAMP TIME	15 min
PRESSURE	9.0 cmH ₂ O

The ramp time value can be set in increments of 5 minutes from 0 (ramp disabled) to the maximum ramp time prescribed, but only when the device is in standby. Press the right control button to raise the setting and the left control button to lower it.

Note: This setting can also be accessed by pressing the ramp button, but only when the device is in standby.

Press the ramp button to access the next parameter on the menu.

SETTINGS	CPAP
COMFORT PRESSURE	6.0 cmH ₂ O
RAMP TIME	15 min
PRESSURE	9.0 cmH ₂ O

The set pressure value does not display unless the Pressure Alter feature is selected. It can be set in increments of 0.5 cmH₂O to the prescribed pressure \pm 1 cmH₂O. Press the right control button to raise the value and the left control button to lower it within the limit from 3 cmH₂O to 20 cmH₂O.

In Bi-level mode

SETTINGS	BILEVEL
COMFORT PRESSURE	6.0 cmH ₂ O
RAMP TIME	15 min

This page displays only when at least one parameter is adjustable. Otherwise, it does not display.

Comfort pressure: Level of pressure produced by the device when the ramp feature starts up.

Ramp time: Set time to reach the inspiratory positive airway pressure (IPAP) and the expiratory positive airway pressure (EPAP) when starting from the comfort pressure.

Note: This page does not display unless the maximum ramp time is other than 0 (see the Parameters screen). The Pressure Alter feature is unavailable in Bilevel mode.

The first parameter on the menu can be accessed by pressing the ramp button. The symbols displayed in the lower part of the screen will change as shown in the following.

SETTINGS	BILEVEL
COMFORT PRESSURE	6.0 cmH ₂ O
RAMP TIME	15 min

The comfort pressure value displays only if the ramp can be activated. The value selected can be set in increments of 0.5 cmH₂O from 3 cmH₂O up to the expiratory positive airway pressure (EPAP), but only when the device is in standby. Press the right control button to raise the setting and the left control button to lower it.

Press the ramp button to access the next parameter on the menu.

SETTINGS	BILEVEL
COMFORT PRESSURE	6.0 cmH ₂ O
RAMP TIME	15 min

The ramp time value can be set in increments of 5 minutes from 0 (ramp disabled) to the maximum ramp time prescribed, but only when the device is in standby. Press the right control button to raise the setting and the left control button to lower it.

Note: This setting can also be accessed by pressing the ramp button, but only when the device is in standby.

Once the settings are made, select the page header by using the ramp button or the on/off button, then press the right control button to access the next page.

3) Parameters

Note: The  symbol indicates that the values displayed on the page cannot be modified.

In CPAP mode

View of the display

PARAMETERS CPAP	
PRESCRIBED P.	8.0 cmH ₂ O
PRESSURE ALTER	<input checked="" type="checkbox"/>
MAX. RAMP TIME	45 min

Prescribed P.: Level of pressure prescribed by your physician.

Pressure Alter: Feature enabled by the physician or home care provider in order to change the prescribed pressure by ± 1 cmH₂O. The symbol indicates that the Pressure Alter feature is enabled, and the symbol indicates that it is disabled.

Max ramp time: Maximum ramp time of from 0 to 45 minutes set by the home care provider.

In Bilevel mode

Views of the display

PARAMETERS BILEVEL	
MAX. RAMP TIME	45 min
EPAP	5.0 cmH ₂ O
IPAP	15.0 cmH ₂ O

Max ramp time: Maximum ramp time of from 0 to 45 minutes set by the home care provider.

EPAP: Set expiratory pressure (low-pressure setting for expiratory positive airway pressure).

IPAP: Set inspiratory pressure (high-pressure setting for inspiratory positive airway pressure).

Press the right control button  to access the next page.

PARAMETERS BILEVEL	
BACKUP FREQUENCY	10 BPM
I/E RATIO	1/2.0
PRESS.RISE TIME	3

Backup frequency: Backup frequency for the device.

I/E ratio: Value of the ratio of inspiratory time to expiratory time.

Press. Rise Time: Pressure rise time set by physician. Time the device takes to rise from the expiratory positive airway pressure (EPAP) to the inspiratory positive airway pressure (IPAP).

Press the right control button  to access the next page.

PARAMETERS BILEVEL	
INSPI.SENSIB.	5
EXPI.SENSIB.	5

Inspi. Sens.: Inspiratory sensitivity set by physician (detection sensitivity of the start of the inspiratory phase of respiration).

Expi. Sens.: Expiratory sensitivity set by physician (detection sensitivity of the start of the expiratory phase of respiration).

Press the right control button  to access the next page.

4) Patient reminders

This screen displays only when the physician or home care provider has set at least one frequency for one type of reminder. Then, the reminder activation date is shown opposite the corresponding reminder.

The display is the same whether the device is set to CPAP mode or Bilevel mode.

Note: The  symbol indicates that the values displayed on the page cannot be modified.

View of the display

REMINDER	
FILTER	09.07.2008
MASK	15.11.2008
SERVICE	09.06.2009
MEMORY CARD	15.08.2008

Filter: Provides a reminder of the date when the filter should be changed.

Mask: Provides a reminder of the date when the mask should be changed.

Service: Provides a reminder of the date when the device pressure should be inspected.

Memory card: Provides a reminder of the date when the memory card should be returned to the home care provider after retrieving the data.

Press the right control button  to access the next page.

5) Setting the clock

The display is the same whether the device is set to CPAP mode or Bilevel mode.

Views of the display

CLOCK	
HOUR	13
MINUTE	45

The time can be set when the device is in standby only.

Press the ramp button to select the hour.

CLOCK	
HOUR	13
MINUTE	45

The value is based on 24 hours and is set in increments of 1. Press the right control button  to raise the setting and the left control button  to lower it.

Press the ramp button to select the minutes.

CLOCK	
HOUR	13
MINUTE	45

The minutes value is based on 60 minutes and is set in increments of 1. Press the right control button  to raise the number of minutes and the left control button  to lower it.

Once the settings are made, select the page header by using the ramp button or the on/off button, then press the right control button  to access the next page.

6) Setting the date

The display is the same whether the device is set to CPAP mode or Bilevel mode.

Views of the display

CALENDAR	
DATE	18
MONTH	9
YEAR	2008

The time can be set when the device is in standby only. Press the ramp button to select the day.

CALENDAR	
DATE	18
MONTH	9
YEAR	2008

Press the right control button  to raise the day setting and the left control button  to lower it.

Press the ramp button to select the month.

CALENDAR	
DATE	18
MONTH	9
YEAR	2008

Press the right control button  to raise the month setting and the left control button  to lower it.

Press the ramp button to select the year.

CALENDAR	
DATE	18
MONTH	9
YEAR	2008


Press the right control button  to raise the year setting and the left control button  to lower it.

Once the settings are made, select the page header by using the ramp button or the on/off button, then press the right control button  to access the next page.

7) Setting the contrast

The display is the same whether the device is set to CPAP mode or Bilevel mode.

Views of the display

Press the ramp button to select the value.

The contrast is set at 5 by default, but can be set from 0 to 10 in increments of 1. Increase the value by using the right control button or decrease it by using the left control button . When the setting is at 10, the contrast is the maximum.

Once your settings are made, select the page header by using the ramp button or the on/off button, then press the on/off button to exit the menu.

General Information

To view the compliance data that has been recorded, press the left control button **(i)** to access this menu.

The data display is identical on the Sandman Duo™ and the Sandman DuoST™.

Note: The  symbol indicates that the values displayed on the page cannot be modified.

Views of the display

GLOBAL INFO	
HOUR COUNTER	86 h 32
TOTAL USE	43 h 25
TOT. SESSION NBR	10
  	

Hour counter: Device operating time.

Total use: Duration of use counter. Total time during which the patient actually breathed in the mask.

Tot. Session Nbr.: Number of start/stop sessions in the device memory.

Press the right control button  to access the next page relating to the last session.

SESSION	08/09/08-16h05
USAGE DURATION	6 h 25
AVERAGE FREQ	15 BPM
AVERAGE LEAK	OK
  	

Usage duration: Length of time during which the patient was actually breathing in the mask.

Average Freq: Average frequency measured during the session.

Average leak: Indicates the average leak, displaying OK, +, ++, +++, or ++++.

The session date and time are displayed in the upper right of the screen.

The pages for the last 5 sessions can be displayed, from the most current to the oldest. Press the right control button  to access the next session.

Press the on/off button to exit the menu.

Reminders about to expire

This feature is the same on the Sandman Duo™ and the Sandman DuoST™.

When one or several reminders are about to expire, the information access symbol in the lower left of the screen changes to the blinking envelope symbol . This appears whether the device is in standby mode (i.e., when the time is displayed) or in operating mode (i.e., when the pressure is displayed).

Pressing the button located opposite the symbol brings up a window with a list of the reminders involved, as the following example shows.

Views of the display

REMINDER	
FILTER	DONE?
MASK	DONE?
SERVICE	DONE?
MEMORY CARD	DONE?
  	

The first parameter is accessible by pressing the ramp button. The symbols displayed in the lower part of the screen change as shown on the screen.

REMINDER	
FILTER	DONE?
MASK	DONE?
SERVICE	DONE?
MEMORY CARD	DONE?
  	

The example asks whether or not the filter has been changed. Press the right control button  to respond yes and the left button  to respond no. If you respond yes, the event reminder disappears from the display.

Once the settings are made, select the page header by using the ramp button or the on/off button, then press the on/off button to exit the menu.

LIST OF AUTHORIZED ACCESSORIES

The following describes the masks, filters and other accessories that may be used with the Sandman Duo™ and Sandman DuoST™ devices.

WARNING

Use authorized accessories only.

Masks

The Sandman Duo™ or Sandman DuoST™ device is designed for use only with nasal interfaces recommended by Puritan Bennett.

Breeze Nasal Interface

P/N: Y-101400-00, Y-101400-XL.

Breeze DreamSeal™ Nasal Interface

P/N: Y-102616-00, Y-103059-00A, Y-103074-00A.

DreamFit Mask (3 sizes)

P/N: DF301, DF302, DF303.

ADAM Interface

P/N: Y-100871-00, Y-100872-00, Y-100873-00, Y-100874-00.

Patient Tubing

Patient Tubing, 22mm diameter, 1.8m (6') length

P/N: M-261000-04.

Patient Tubing, 22mm diameter, 2.4m (8') length

P/N: S-776274-00C.

Other Accessories

Oxygen Adapter

P/N: Y-616155-03B.

Battery Cable

P/N: M-214530-00

Cigarette-lighter Cable

P/N: M-214530-01

RS232 cable, 2 m Sandman CPAP/Bilevels

P/N: M-214831-05

RS232 cable, 15 m Sandman CPAP/Bilevels

P/N: M-214831-04

USB cable, 2 m Sandman CPAP/Bilevels

P/N: M-214831-06

Memory Card

P/N: M-214830-00

Power Supply Cord

P/N: M-660601-43

Sandman CPAP/Bilevels Cover

P/N: M-314830-06

Air Intake Filter (foam, reusable)

P/N: M-314850-08 (x10); M-414840-06 (x96)

Optional High-efficiency Filter (disposable)

P/N: M-314850-09 (x10); M-414841-07 (x100)

Filter holder kit

P/N: M-314820-02

Water Chamber Kit (includes base, cover and gasket)

P/N: M-314810-15

Memory card reader

P/N: M-214830-03

Carrying Case

P/N: M-814805-00

KnightControl remote control

P/N: M-114700-US

CLEANING AND MAINTENANCE

To ensure reliable and efficient use of the **Sandman Duo™** or **Sandman DuoST™** devices, follow the maintenance instructions shown below.

Sandman Duo™ or Sandman DuoST™ devices

The device casing does not require any special maintenance. If necessary, clean the device casing by using a cloth moistened with warm soapy water (for example, mild liquid soap). Dry the device thoroughly.

WARNINGS

- Always unplug the device from the electrical source before cleaning. Let it cool down.
- Do not use sprays. Chemical product residue could enter the air outlet, the air intake filter, or the device interior, causing airway irritation.
- Do not immerse the device in water.
- Never allow water to spill onto the device or to drip into any of its openings.
- Never use the device without making sure that an air intake filter is installed.
- Do not insert foreign objects into the device.

CAUTION

The Sandman Duo™ or Sandman DuoST™ device should be opened by skilled personnel only.

At least once a year, contact your home care provider to have the device pressure inspected.

Air intake filters

Change the filters periodically, when they are torn or soiled (see Item 8 of Figure 5 on page 7) or when the device displays a reminder to do so.

Washable filter

Clean the reusable air-intake filter at least once a week, and more often if necessary, as follows:

1. Wash the filter in hot soapy water.
2. Rinse well to eliminate any trace of soap.
3. Press dry in a terry towel.
4. Allow to dry completely before reinstalling. Do not use a filter that is not completely dry.
5. Reinstall the filter at the rear of the device.

Optional filter

The high-efficiency filter, which is offered as an option, cannot be washed and therefore must be discarded. It must be changed at least once a month and more frequently if visibly dirty.

Cleaning the water chamber

1. Remove the Sandman Duo™ or Sandman DuoST™ patient circuit. Turn the device off and allow it to cool.
2. To disassemble the water chamber, see "" on page 44.
3. Remove any remaining water.
4. Clean the silicone gasket and both parts of the water chamber daily in warm soapy water and rinse well with tap water. The two parts of the water chamber without the silicone gasket may also be washed in a dishwasher. Allow to dry and wipe the exterior with a clean cloth.
5. Once a week, soak the water chamber for 15 minutes in a solution of 1 part white vinegar to 9 parts water. Rinse well with tap water. Dry the exterior with a clean cloth.
6. To prevent bacterial growth, do not allow water in the reservoir to stagnate.

WARNING

Check that the heating element is completely dry before reconnecting the device to the power supply.

Patient tubing

Follow the instructions in the directions for use that came with the patient tubing.

Mask

Follow the instructions in the directions for use that came with the mask.

Oxygen adapter

Follow the instructions in the directions for use that came with the oxygen adapter.

TROUBLESHOOTING




Helpful hints

Problem	Possible cause	Suggestion
Your nose is cold.	The room temperature is too low.	Raise the room temperature.
	The delivered air is too cold.	Place patient tubing under a cover and run next to your body to reduce heat loss.
Runny Nose.	Reaction to air flow and pressure levels	Contact your medical technical team or your attending physician.
Nose or throat is dry or irritated.	The air is too dry.	Use the integrated humidifier. Raise the humidity level following the instructions under "", on page 44.
	There is no water in the water chamber.	Check the water level in the water chamber. If necessary, raise the water level following the instructions under "", on page 44.
Nose, sinus or ear pain.	Sinus infection or nasal congestion.	Contact the attending physician immediately.
Redness of skin in contact with the mask.	Possible pressure sores indicating the headgear is too tight or the size is too small.	Adjust the headgear. Contact the attending physician or home care provider to try another size mask.
	Allergic reaction to mask components.	Stop using the mask. Contact the attending physician or home care provider.
Dry or irritated eyes.	Air leakage around the mask.	Reposition the mask. Ask the attending physician or home care provider to try another size mask.
The air delivered by the device is too hot.	Air intake filters may be dirty. Air intake is clogged.	Clean or replace filters as necessary (see "Cleaning and Maintenance", page 32). Move away from all linens and clothing.
	Room temperature is too high.	Lower the room temperature. Make sure that the Sandman Duo or Sandman DuoST device is well away from any source of heat. Remove patient tubing from cover.
Discomfort from feeling that the pressure is too high.	Device pressure.	Getting use to the nasal pressure takes some time. Use the pressure ramp when you go to sleep (see "", page 44). Relax and breathe slowly through the nose.
	The Sandman Duo™ or Duo™ST device is set to Bilevel mode.	If the pressure delivered by the device seems to have changed, contact the home care provider to have the device pressure checked.
The device does not display the correct pressure.	The ramp is enabled.	The pressure levels have been prescribed by your physician and can be changed only by medical prescription.
		Confirm that the ramp indicator is displayed. Disable the ramp feature following the procedure described under "", page 44.

Problem	Possible cause	Suggestion
Symptoms of sleep apnea syndrome recur.	The device is not set to the correct pressure or is not operating properly. Your physical condition or pressure needs may have changed.	Ask your home care provider to check the device operation. Contact your attending physician.
The Sandman Duo™ or Sandman DuoST™ device does not light up (no display).	The power supply cord is not correctly plugged into the power port. No power. The device's internal fuse is defective.	Check the connections between the device, the power cord, and the power source. Plug another device (e.g., lamp, radio, etc.) into the power source to confirm that current is available at the source. Contact your home care provider.
The Sandman Duo™ or Sandman DuoST™ device is not operating smoothly and appears to be intermittent.	Very high electromagnetic interference.	Move the device away from sources of interference, such as halogen lamps, cordless phones, etc.
Humidifier does not seem to heat water.	Humidification is set at 0. The heating element is defective.	When the device is battery-powered, the heated humidification feature cannot be used. Check that the device is equipped with an integrated heated humidifier, and humidification is set at a level from 1 to 10 according to the instructions under "", on page 44. Confirm that the heated humidifier indicator is displayed on the screen and is not blinking. If the indicator is blinking or is not present, contact your home care provider.
Water spilled into the device.		Unplug the device and allow it to dry at least 24 hours. Re-apply power and then check for proper operation.
Excessive droplets of water appear in the patient circuit or in the mask.	The water chamber water level is too high. A small amount of water droplets are normal, especially during winter months. Excessive condensation of water vapor.	Check that the water level is not above the water chamber maximum fill level; if it is, pour out the excess water. Lower the humidity level following the instructions under "", on page 44. Route the patient circuit under your blanket.

Problem	Possible cause	Suggestion
White, powdery deposits appear in the water chamber.	Traces of scale from the tap water.	Rub the traces with a sponge and liquid dish detergent. Soak the water chamber in a solution of 1 part white vinegar to 9 parts water as described under "Cleaning the water chamber", on page 33. Rinse well with tap water. Dry the exterior with a clean cloth.

Device messages

Code	Description	Action corrective
MASK UNPLUGGED is displayed on the screen.	The mask is disconnected.	Check the connections between the mask, the patient circuit, and the device. This message disappears when any of the four buttons on the front of the device is pressed.
HIGH FREQUENCY is displayed on the screen.	An overly high frequency has been detected.	Press the device ON/OFF button twice. Check the patient circuit connections. If the problem recurs, contact your home care provider.
 is blinking in the center of the screen.	The device has detected an operating error.	Unplug the device from all sources of electrical power. Plug it in again and turn on. If the error persists, contact your home care provider.
 The symbol  is blinking in the center of the humidification level settings page.	The device has detected a humidifier operating error.	If the device is connected to a battery, the heated humidification feature cannot be used and the symbol blinks. Unplug the device from all sources of electrical power. Plug it in again and turn on. If the error persists, contact your home care provider.

TECHNICAL CHARACTERISTICS

Sandman Duo™ or Sandman DuoST™ Device Performance

Device pressure range:	3 cmH ₂ O to 20 cmH ₂ O in CPAP mode 3 cmH ₂ O to 25 cmH ₂ O in Bilevel mode Device adjustable in increments of 0.5 cmH ₂ O
Maximum achievable pressure at the patient connection port under single fault condition:	30 cmH ₂ O in CPAP mode or in Bilevel mode if IPAP = EPAP 40 cmH ₂ O in Bilevel mode if IPAP ≠ EPAP
Maximum pressure that can be set:	20 cmH ₂ O in CPAP mode 25 cmH ₂ O in Bilevel mode
Ramp duration:	0 to 45 minutes ± 1 minute Device adjustable in increments of 5 minutes.
Sound level measured in accordance with ISO 17510-1:2002:	< 26 dBA
Absolute maximum sound level at one meter:	30 dBA at 10 cmH ₂ O without water chamber
Battery life of mainboard:	5 years
Patient-side connection aperture:	Tubing of 22mm external diameter compliant with the EN1281 standard
Air intake filter efficiency:	Optional high-efficiency filter, disposable fabric HEPA paper, which is 90% effective for particles < 3 microns. High Density Foam (reusable) and is 90% effective with particles > 20 microns.

Integrated Humidifier Performance on Sandman Duo™ HC and Sandman DuoST™ HC

Humidification rate:	10 mg H ₂ O/s minimum at maximum setting
Heating time:	≤ 45 minutes
Pressure loss as a function of flow:	0.5 cmH ₂ O at 1 L/sec
Maximum service pressure:	25 cmH ₂ O
Maximum gas temperature delivered from patient tubing:	41°C (+106 °F)

Physical characteristics

Dimensions (W x H x D):	6.9 x 4.7 x 7.5 in approximately
Carrying case Dimensions (W x H x D):	11.0 x 5.5 x 13.8 in approximately
Weight:	2.2 lbs approximately without humidifier 2.64 lbs approximately with humidifier
Operating volume: (minimal gas volume in water chamber)	230 mL
Useful volume of water chamber:	350 mL
Maximum liquid fill level:	Indicated on water chamber

Transportation and storage conditions

Relative pressure range:	500 hPa to 1060 hPa (1 hecto-Pascal ~ 1 cmH ₂ O)
Temperature range:	-20°C to +60°C (-4 °F to +140 °F)
Relative humidity:	Up to 95% without condensation

Operating conditions for the device and humidifier

Pressure range:	730 hPa to 1060 hPa (1 hecto-Pascal ~ 1 cmH ₂ O)
Temperature:	+5°C to +40°C (+41 °F to +104 °F) without the humidifier
Temperature:	+5°C to +35°C (+41 °F to +95 °F) with the humidifier
Relative humidity:	between 10% and 95% without condensation
Altitude range:	approximately 0–9,000 ft

Device electrical characteristics

Input voltage (auto-ranging):	100, 115, 120, 230, and 240 VAC at 50 or 60Hz
Universal internal electrical power with the detachable electrical power cord (which varies by country).	
Maximum power consumption:	90VA at maximum setting.
Power consumption is based on the type of model. Maximum consumption is to be considered as 90VA.	
Model Sandman Duo™ or Sandman DuoST™ –HC is considered as the highest consumption unit.	
Current consumed at 20cmH ₂ O with a 4mm leak:	130 mA eff. at 230 VAC Direct current

Meaning of symbols displayed on the screen

Symbol	Description	Symbol	Description
	Integrated humidifier present, heater disabled.		Heated humidification. The symbol blinks when the heated humidifier has a problem or when the device is connected to a battery.
(i)	Access to device information.		Access to device settings.
	Parameters displayed on the screen can be changed.		Parameters displayed on the screen are for information and cannot be changed.
	Feature enabled.		Feature disabled.
	YES.		NO.
	Provides access to the next page available for display on the screen.		Provides access to the preceding page displayed on the screen.
	Allows the value of the parameter displayed on the screen to be raised.		Allows the value of the parameter displayed on the screen to be lowered.
	The symbol blinks when the ramp feature is selected but still disabled.		Indicates that compliance is being recorded.
C	CPAP operating mode.		Bilevel operating mode (two pressure levels).
	The symbol blinks when one or several patient reminders are about to expire.		Memory card present. The symbol blinks if the memory card is not operational.
	Indicates that a remote control is connected.		Indicates that a battery is connected.
	The device has detected an operating error.		

Symbols on the device

On/off button.

Access ramp time setting.

Raise setting value

Lower setting value

IPX1

Device protected against falling drops of water.

In compliance with the requirements of European Directive 2002/96/EC on Waste Electrical and Electronic Equipment (WEEE), used devices are to be collected separately from waste. Contact your home care provider.

Class II device.

Type B device.


Consult use instructions before using this device.

Rx ONLY

For the United States only:


Federal law restricts this device to sale by or on the order of a physician.

On packaging: this icon means "This side up", showing how to place the package before opening.

On packaging: this icon means "Keep dry", as the package needs to be protected against moisture and water

On packaging: this icon means "Fragile", as the package needs to be handled with care

Standards Compliance

The Sandman Duo™ or Sandman DuoST™ device complies with the requirements of the following standards:

- IEC 60 601-1:1988 + A1:1991 + A2:1995 + US/CAN deviations: Medical electrical equipment. Part 1: General requirements for safety.
- IEC 60 601-1-2:2001 + A1: 2004: Medical electrical equipment. Part 1: General requirements for safety. Collateral standard: Electromagnetic compatibility – Requirements and tests.
- EN ISO 17510-1:2002: Sleep Apnoea Breathing Therapy Devices.
- EN ISO 8185: 2007: General Requirements for Humidification Systems.
- European Parliament Directive 2002/95/EC on the restriction of the use of certain hazardous substances (RoHS) in electrical and electronic equipment.

WARNING

Devices connected via a serial RS232/USB connection must be EN60950 compliant.

Essential performance requirements for electromagnetic compatibility

Since the Sandman Duo™ and Sandman DuoST™ are not life support systems, they are vulnerable to disruption during EMC disturbances. However, these devices should resume normal operations once the disturbance is past.

Electromagnetic emissions			
<i>Sandman Duo™ / Sandman DuoST™</i> devices are intended for use in the electromagnetic environment specified below. The customer or the user of a <i>Sandman Duo™ / Sandman DuoST™</i> device should assure that it is used in such an environment.			
Emission test	Conformity	Recommended electromagnetic environment	
RF emissions CISPR 11	Group 1	<i>Sandman Duo™ / Sandman DuoST™</i> devices uses RF energy only for its internal operations. Therefore its RF emissions are very low and cannot be assumed to interfere with nearby electronic equipment.	
RF emissions CISPR 11	Class B	<i>Sandman Duo™ / Sandman DuoST™</i> devices are suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.	
Harmonic emissions IEC 61000-3-2	Class A		
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies		
Electromagnetic Immunity			
<i>Sandman Duo™ / Sandman DuoST™</i> devices are intended for use in the electromagnetic environment specified below. The customer or the user of a <i>Sandman Duo™ / Sandman DuoST™</i> device should assure that it is used in such an environment.			
Immunity test	Test level IEC 60601	Compliance level	Recommended electro-magnetic environment
Electro-static discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input / output lines	± 2 kV for power supply lines ± 1 kV for input / output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV lines / lines ± 2 kV lines/earth	± 1 kV lines / lines ± 2 kV lines/earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5% U _T (>95% dip in U _T) for 0.5 cycles 40% U _T (60% dip in U _T) for 5 cycles 70% U _T (30% dip in U _T) for 25 cycles	< 5% U _T (>95% dip in U _T) for 0.5 cycles 40% U _T (60% dip in U _T) for 5 cycles 70% U _T (30% dip in U _T) for 25 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the <i>Sandman Duo™ / Sandman DuoST™</i> device requires continued operation during power mains, it is recommended that the <i>Sandman Duo™ / Sandman DuoST™</i> device be powered from an uninterruptible power supply or a battery.

Electromagnetic Immunity (continued)°

	< 5% U_T (>95% dip in U_T) for 5 s	< 5% U_T (>95% dip in U_T) for 5 s	
Magnetic fields in the Power Supply Frequencies (50/60 Hz) IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note: U_T is the AC mains voltage prior to application of the test level

Electromagnetic Immunity - conducted and radiated RF

Sandman Duo™ / Sandman DuoST™ devices are intended for use in the electromagnetic environment specified below. The customer or the user of the **Sandman Duo™ / Sandman DuoST™** device should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms 150 kHz to 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the Sandman Duo™ / Sandman DuoST™ device, including cables. Recommended separation distance : $d = 1.2\sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2\sqrt{P} \text{ 80 MHz to 800 MHz}$ $d = 2.3\sqrt{P} \text{ 800 MHz to 2.5 GHz}$

Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).

Fields strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b.

Interference may occur in the vicinity of equipment marked with the following symbol:

NOTE 1: At 80 MHz and 800 MHz, the highest frequency range applies

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

^{a)} Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which a **Sandman Duo™ / Sandman DuoST™** device is used exceeds the applicable RF compliance level above, the **Sandman Duo™ / Sandman DuoST™** device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the **Sandman Duo™ / Sandman DuoST™** device.

^{b)} Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V / m.

Recommended separation distance between portable and mobile RF communications equipment and a *Sandman Duo™ / Sandman DuoST™* device

Sandman Duo™ / Sandman DuoST™ devices are designed to be used in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of a *Sandman Duo™ / Sandman DuoST™* device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the *Sandman Duo™ / Sandman DuoST™* device as recommended below, according to the maximum power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance of the highest frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

Compliant cables

Warning

The use of accessories and cables other than those specified, with the exception of parts sold by Puritan Bennett as replacements for internal components, may result in increased emissions or decreased immunity of the *Sandman Duo™ / Sandman DuoST™* devices.

M-214831-05 SANDMAN CPAP RS232 2m cable	6.5 ft (2m)
M-214831-04 SANDMAN CPAP RS232 15m cable	50 ft (15m)
M-214831-06 SANDMAN CPAP USB 2m cable	6.5 ft (2m)

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