

HOW TO IDENTIFY HUMIDIFIER CHAMBER

Humidification is an important part of your sleep therapy treatment to aid in comfort of use. Your Respirationics M Series sleep therapy system includes a removable water chamber that should be replaced as instructed by your Homecare Provider. Since the release of the M Series sleep system, Respirationics has enhanced the replaceable water chamber unit to enable easier use and a better seal. The following pictures will help you identify for your homecare provider one of three water chamber models that is compatible with your humidification platform.

With your current chamber in hand, please order from your Homecare Provider, the part number that corresponds with the pictures in one of the three columns.



BLUE LATCH



BLACK LATCH



BLACK LATCH



BLUE MIDDLE SEAL




BLACK MIDDLE SEAL



BLACK MIDDLE SEAL



**SIDE-BY-SIDE
OPENINGS**



**SIDE-BY-SIDE
OPENINGS**



**STAGGERED
OPENINGS**

R2B Chamber Kit
Part Number 1003757

R2B Chamber Kit
Part Number 1003757

R3B Chamber Kit
Part Number 1003758