

Mask settings for System One Resistance Control

System One Resistance Control is a unique technology incorporated into the Philips Respironics Sleep Therapy System that knows how to optimize the device's pressure delivery for any Philips Respironics mask. The table below outlines the System One Resistance Control settings for Philips Respironics' nasal, minimal contact, and full-face masks.

Nasal masks	Size	System One setting
ComfortGel	All sizes	X1
ComfortFusion	All sizes	X1
ComfortSelect	All sizes	X1
ComfortClassic	All sizes	X1
Profile Lite	All sizes	X1
Minimal contact masks	Size	System One setting
OptiLife	Pillows – petite	X3
	Pillows – small	X4
	Pillows – medium	X2
	Pillows – large	X1
	CradleCushion all sizes	X1
ComfortLite 2	Direct Seal – 1	X2
	Direct Seal – 2	X2
	Direct Seal – 3	X2
	Direct Seal – 4	X1
	Direct Seal – 5	X2
	Direct Seal – 6	X2
	Pillows – petite	X4
	Pillows – small	X4
	Pillows – medium	X2
	Pillows – large	X1
Simple Cushion all sizes	X1	
Full-face masks	Size	System One setting
FitLife	All sizes	X1
FullLife	All sizes	X1
ComfortGel Full	All sizes	X1
ComfortFull 2	All sizes	X1

PHILIPS

RESPIRONICS

**Philips Healthcare is part of
Royal Philips Electronics**

How to reach us

www.philips.com/healthcare

healthcare@philips.com

fax: +31 40 27 64 887

Asia

+852 2821 5888

Europe, Middle East, Africa

+49 7031 463 2254

Latin America

+55 11 2125 0744

North America

+1 425 487 7000

800 285 5585 (toll free, US only)

Philips Respironics
1010 Murry Ridge Lane
Murrysville, PA 15668

Customer Service

+1 724 387 4000

800 345 6443

(toll free, US only)

Respironics Asia Pacific

+65 6298 1088

Respironics Australia

+61 (2) 9666 4444

Respironics Deutschland

+49 8152 93 06 0

**Respironics Europe,
Middle East, Africa**

+33 1 47 52 30 00

Respironics France

+33 2 51 89 36 00

Respironics Italy

+39 03 62 63 43 1

Respironics Sweden

+46 8 120 45 900

Respironics Switzerland

+41 6 27 45 17 50

Respironics United Kingdom

+44 800 1300 845

www.philips.com/respironics

Respironics, ComfortClassic, ComfortFull, ComfortFusion, ComfortGel, ComfortLite, ComfortSelect, FitLife, FullLife, OptiLife, and Profile Lite are trademarks of Respironics, Inc. and its affiliates. All rights reserved.



©2009 Koninklijke Philips Electronics N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

CAUTION: US federal law restricts these devices to sale by or on the order of a physician.

Hoech JJ 09/24/09 MCI 4102845 PN 1067419